



# SAFETY DATA SHEET

Revision date 01-Jan-2021  
Version 1.0

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

### Product identifier

Product Name Acid green 50  
Product Code S10108

### Recommended use of the chemical and restrictions on use

For research use only. Not intended for diagnostic or therapeutic use.

### Details of the supplier of the safety data sheet

DAWN SCIENTIFIC  
121 Liberty Street, Metuchen, NJ  
1-800-DAWN-SCI | 732-902-6300

### Emergency telephone number

Chemtrec  
1.800.424.9300 (Within USA)  
+1.703.527.3887 (Outside USA)

## 2. HAZARDS IDENTIFICATION

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

### Classification

Acute toxicity - Oral	Category 4
Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2A
Specific target organ toxicity (single exposure)	Category 3

### Label elements

Signal word	Warning
Hazard statements	Harmful if swallowed Causes skin irritation Causes serious eye irritation May cause respiratory irritation

Symbols/Pictograms



Precautionary Statements - Prevention

Wash face, hands and any exposed skin thoroughly after handling  
Do not eat, drink or smoke when using this product  
Wear protective gloves/protective clothing/eye protection/face protection  
Avoid breathing dust/fume/gas/mist/vapors/spray  
Use only outdoors or in a well-ventilated area

Precautionary Statements - Response

IF exposed or concerned: Get medical advice/attention  
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing  
If eye irritation persists: Get medical advice/attention  
IF ON SKIN: Wash with plenty of soap and water  
If skin irritation occurs: Get medical advice/attention  
Take off contaminated clothing and wash before reuse  
IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing  
Call a POISON CENTER or doctor/physician if you feel unwell  
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell  
Rinse mouth

Precautionary Statements - Storage

Store in a well-ventilated place. Keep container tightly closed  
Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

**Hazards not otherwise classified (HNOC)**

Hazards not otherwise classified (HNOC)

Not applicable

**Other Information**

**NFPA**

Health hazards 2  
Flammability 0  
Stability 0  
Physical and chemical properties -



**HMIS**

Health hazards 2  
Flammability 0  
Physical hazards 0  
Personal protection -

**3. COMPOSITION/INFORMATION ON INGREDIENTS**

CAS No 3087-16-9  
Molecular Weight 576.62  
Formula C<sub>27</sub>H<sub>25</sub>N<sub>2</sub>O<sub>7</sub>S<sub>2</sub>•Na

Chemical Name	CAS No	Weight %	Oral LD50	Dermal LD50	Inhalation LC50
Acid Green 50	3087-16-9	>98	= 2000 mg/kg ( Rat )	-	-

**4. FIRST AID MEASURES**

**First Aid Measures**

General advice If symptoms persist, call a physician. Do not breathe dust/fume/gas/mist/vapors/spray. Do not get in eyes, on skin, or on clothing.

Eye contact Immediately flush with plenty of water. After initial flushing, remove any contact lenses and continue flushing for at least 15 minutes. Keep eye wide open while rinsing. If symptoms persist, call a physician.

Skin Contact Consult a physician if necessary. Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. Wash contaminated clothing before reuse. Wash off immediately with plenty of water. If skin irritation persists, call a physician. Immediate medical attention is not required.

Inhalation Remove to fresh air Call a physician If breathing is irregular or stopped, administer artificial respiration Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation Artificial respiration and/or oxygen may be necessary Immediate medical attention is not required Move to fresh air in case of accidental inhalation of vapors If symptoms persist, call a physician

---

Ingestion	Rinse mouth. Drink plenty of water. If symptoms persist, call a physician. Do NOT induce vomiting. Clean mouth with water and drink afterwards plenty of water. Never give anything by mouth to an unconscious person. Call a physician.
Self-protection of the first aider	Use personal protective equipment as required.

**Most important symptoms and effects, both acute and delayed**

Symptoms No information available.

**Indication of any immediate medical attention and special treatment needed**

Note to physicians Treat symptomatically.

---

**5. FIRE-FIGHTING MEASURES**

**Suitable Extinguishing Media**

Suitable Extinguishing Media Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media None.

**Specific hazards arising from the chemical**

Specific hazards arising from the chemical No information available.

Hazardous combustion products Hydrogen sulfide. Carbon oxides. Nitrogen oxides (NOx).

**Explosion data**

Sensitivity to Mechanical Impact No information available.

Sensitivity to Static Discharge No information available.

**Protective equipment and precautions for firefighters**

Protective equipment and precautions for firefighters As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

---

**6. ACCIDENTAL RELEASE MEASURES**

**Personal precautions, protective equipment and emergency procedures**

Personal precautions Use personal protective equipment as required. Evacuate personnel to safe areas. Keep people away from and upwind of spill/leak.

**Environmental precautions**

Environmental precautions Prevent entry into waterways, sewers, basements or confined areas. Do not flush into surface water or sanitary sewer system. Prevent further leakage or spillage if safe to do so. Prevent product from entering drains. See Section 12 for additional Ecological Information.

**Methods and material for containment and cleaning up**

Methods for containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Use personal protective equipment as required. Cover powder spill with plastic sheet or tarp to minimize spreading and keep powder dry. Take up mechanically, placing in appropriate containers for disposal. Avoid creating dust. Clean contaminated surface thoroughly. Soak up with inert absorbent material. Dam up. Pick up and transfer to properly labeled containers. Sweep up and shovel into suitable containers for disposal. After cleaning, flush away traces with water. Take precautionary measures against static discharges.

---

**7. HANDLING AND STORAGE**

---

### Precautions for safe handling

Advice on safe handling

Avoid contact with skin, eyes or clothing. Use personal protective equipment as required. Wash contaminated clothing before reuse. Do not breathe dust/fume/gas/mist/vapors/spray. Do not eat, drink or smoke when using this product. Use with local exhaust ventilation. Handle in accordance with good industrial hygiene and safety practice. Corrosive hazard. Wear protective gloves/clothing and eye/face protection.

### Conditions for safe storage, including any incompatibilities

Storage Conditions

Keep container tightly closed in a dry and well-ventilated place. Keep out of the reach of children. Keep containers tightly closed in a cool, well-ventilated place. Keep in properly labeled containers. Store at room temperature.

Incompatible materials

None known based on information supplied.

---

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### Control parameters

Exposure Guidelines

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

### Appropriate engineering controls

Engineering Controls

Showers  
Eyewash stations  
Ventilation systems

### Individual protection measures, such as personal protective equipment

Eye/face protection

Tight sealing safety goggles. Face protection shield.

Skin and Body Protection

Wear protective gloves and protective clothing.

Respiratory protection

If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

General Hygiene Considerations

When using do not eat, drink or smoke. Regular cleaning of equipment, work area and clothing is recommended.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

Physical State

Solid

Appearance

powder

Odor

No information available

### Property

### Values

pH

No information available

Melting point/freezing point

No information available

Boiling point

No information available

Flash point

No information available

Density

No information available

Evaporation rate

No information available

Upper flammability limits

No information available

Lower flammability limit

No information available

Vapor pressure

No information available

Vapor density

No information available

Specific gravity

No information available

Water solubility

No information available

Solubility in other solvents

No information available

Partition coefficient

No information available

---

Autoignition temperature	No information available
Decomposition temperature	No information available
Kinematic viscosity	No information available
Explosive properties	No information available
Oxidizing properties	No information available

## 10. STABILITY AND REACTIVITY

---

Reactivity	Not applicable
Chemical stability	Stable under recommended storage conditions.
Possibility of Hazardous Reactions	None under normal processing.
Hazardous polymerization	No information available.
Conditions to avoid	Extremes of temperature and direct sunlight.
Incompatible materials	Strong oxidizing agents.
Hazardous Decomposition Products	Hydrogen sulfide. Carbon oxides. Nitrogen oxides (NOx).

## 11. TOXICOLOGICAL INFORMATION

---

### Information on likely routes of exposure

Inhalation	No data available.
Eye contact	No data available.
Skin Contact	No data available.
Ingestion	No data available.

### Information on toxicological effects

Symptoms	No information available.
----------	---------------------------

### Delayed and immediate effects as well as chronic effects from short and long-term exposure

Chronic Toxicity	Avoid repeated exposure.
------------------	--------------------------

### Numerical measures of toxicity - Product Information

Unknown acute toxicity	No information available
The following values are calculated based on chapter 3.1 of the GHS document	
ATEmix (oral)	2000 mg/kg

## 12. ECOLOGICAL INFORMATION

---

Ecotoxicity	May cause long lasting harmful effects to aquatic life
-------------	--

100% of the mixture consists of components(s) of unknown hazards to the aquatic environment.

Persistence and degradability	No information available.
Bioaccumulation	No information available.
Mobility	No information available.

## 13. DISPOSAL CONSIDERATIONS

---

Disposal of wastes	Disposal should be in accordance with applicable regional, national and local laws and regulations.
Contaminated packaging	Do not reuse container.

---

Other Information

Waste codes should be assigned by the user based on the application for which the product was used.

## 14. TRANSPORT INFORMATION

---

DOT Not regulated

IMDG Not regulated

IATA Not regulated

## 15. REGULATORY INFORMATION

---

### International Inventories

All of the components in the product are on the following inventory lists

TSCA (United States): Canada (DSL/NDSL) Europe (EINECS/ELINCS/NLP) Australia (AICS) China (IECSC) ENCS (Japan):

Chemical Name	TSCA	DSL	NDSL	EINECS	ELINCS	ENCS	IECSC	KECL	PICCS	AICS
Acid Green 50	X	X	-	X	-	X	X	-	-	X

X - Listed

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

### US Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

#### **SARA 311/312 Hazard Categories**

Acute health hazard	Yes
Chronic Health Hazard	No
Fire hazard	No
Sudden release of pressure hazard	No
Reactive hazard	No

#### **CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

### US State Regulations

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

#### **U.S. State Right-to-Know Regulations**

This product does not contain any substances regulated by state right-to-know regulations

---

## 16. OTHER INFORMATION

---

Revision note

No information available

### **Disclaimer**

The information provided in this Material Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.