



# SAFETY DATA SHEET

Revision date 01-Jan-2021  
Version 1.0

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

### Product identifier

Product Name Nile Red  
Product Code S10134

### Recommended use of the chemical and restrictions on use

For research use only. Not intended for diagnostic or therapeutic use.

### Details of the supplier of the safety data sheet

DAWN SCIENTIFIC  
121 Liberty Street, Metuchen, NJ  
1-800-DAWN-SCI | 732-902-6300

### Emergency telephone number

Chemtrec  
1.800.424.9300 (Within USA)  
+1.703.527.3887 (Outside USA)

## 2. HAZARDS IDENTIFICATION

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.122).

### Classification

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS)

### Label elements

Signal word Not classified  
Hazard statements Not classified  
Symbols/Pictograms Not classified

Precautionary Statements - Prevention  
Precautionary Statements - Response

Wash hands thoroughly after handling  
IF exposed or concerned: Get medical advice/attention

### Hazards not otherwise classified (HNOC)

Hazards not otherwise classified (HNOC) Not applicable

### Other Information

Unknown acute toxicity 100% of the mixture consists of ingredient(s) of unknown toxicity.

### NFPA

Health hazards 0  
Flammability 1  
Stability 0  
Physical and chemical properties -



### HMIS

Health hazards 0  
Flammability 1  
Physical hazards 0  
Personal protection -

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

CAS No 7385-67-3  
Molecular Weight 318.37  
Formula  $C_{20}H_{18}N_2O_2$

Chemical name	CAS No.	Weight-%	Oral LD50	Dermal LD50	Inhalation LC50
---------------	---------	----------	-----------	-------------	-----------------

Nile Red	7385-67-3	>98	-	-	-
----------	-----------	-----	---	---	---

## 4. FIRST AID MEASURES

### First Aid Measures

General advice	Consult a physician if necessary. Remove to fresh air.
Eye contact	Wash with plenty of water.
Skin Contact	Wash skin with soap and water.
Inhalation	Remove to fresh air If breathing is difficult, give oxygen If not breathing, give artificial respiration
Ingestion	Never give anything by mouth to an unconscious person. Clean mouth with water.

### Most important symptoms and effects, both acute and delayed

Symptoms	No information available.
----------	---------------------------

### Indication of any immediate medical attention and special treatment needed

Note to physicians	Treat symptomatically.
--------------------	------------------------

## 5. FIRE-FIGHTING MEASURES

### Suitable Extinguishing Media

Suitable Extinguishing Media	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
Unsuitable Extinguishing Media	No information available.

### Specific hazards arising from the chemical

Specific hazards arising from the chemical	No information available.
Hazardous combustion products	Carbon oxides. Nitrogen oxides (NOx).

### Explosion data

Sensitivity to Mechanical Impact	No information available.
Sensitivity to Static Discharge	No information available.

### Protective equipment and precautions for firefighters

Protective equipment and precautions for firefighters	As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.
---	--

## 6. ACCIDENTAL RELEASE MEASURES

### Personal precautions, protective equipment and emergency procedures

Personal precautions	Ensure adequate ventilation, especially in confined areas.
----------------------	--

### Environmental precautions

Environmental precautions	See Section 12 for additional Ecological Information.
---------------------------	---

### Methods and material for containment and cleaning up

Methods for containment	Prevent further leakage or spillage if safe to do so.
Methods for cleaning up	Use personal protective equipment as required. Take up mechanically, placing in appropriate containers for disposal. Clean contaminated surface thoroughly. Cover powder spill with plastic sheet or tarp to minimize spreading and keep powder dry. Avoid creating dust.

---

## 7. HANDLING AND STORAGE

---

### Precautions for safe handling

Advice on safe handling Corrosive hazard. Wear protective gloves/clothing and eye/face protection. Handle in accordance with good industrial hygiene and safety practice.

### Conditions for safe storage, including any incompatibilities

Storage Conditions Keep containers tightly closed in a dry, cool and well-ventilated place. Store at room temperature.  
Incompatible materials None known based on information supplied.

---

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

---

### Control parameters

Exposure Guidelines This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

### Appropriate engineering controls

Engineering Controls Showers  
Eyewash stations  
Ventilation systems

### Individual protection measures, such as personal protective equipment

Eye/face protection Wear safety glasses with side shields (or goggles).  
Skin and Body Protection Wear protective gloves and protective clothing.  
Respiratory protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.  
General Hygiene Considerations Handle in accordance with good industrial hygiene and safety practice.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

---

Physical State Solid  
Appearance powder  
Odor No information available

### Property

<u>Property</u>	<u>Values</u>
pH	No information available
Melting point/freezing point	203 °C
Boiling point	484.7 °C
Flash point	246.9 °C
Liquid Density	1.2 g/cm <sup>3</sup>
Evaporation rate	No information available
Upper flammability limits	No information available
Lower flammability limit	No information available
Vapor pressure	0.0 mmHg
Vapor density	No information available
Specific gravity	No information available
Water solubility	No information available
Solubility in other solvents	No information available
Partition coefficient	4.38
Autoignition temperature	No information available
Decomposition temperature	No information available

---

Kinematic viscosity	No information available
Explosive properties	No information available
Oxidizing properties	No information available

## 10. STABILITY AND REACTIVITY

---

Reactivity	Not applicable
Chemical stability	Stable under recommended storage conditions.
Possibility of Hazardous Reactions	None under normal processing.
Hazardous polymerization	No information available.
Conditions to avoid	Extremes of temperature and direct sunlight.
Incompatible materials	Strong oxidizing agents.
Hazardous Decomposition Products	Carbon oxides. Nitrogen oxides (NOx).

## 11. TOXICOLOGICAL INFORMATION

---

### Information on likely routes of exposure

Inhalation	No data available.
Eye contact	No data available.
Skin Contact	No data available.
Ingestion	No data available.

### Information on toxicological effects

Symptoms	No information available.
----------	---------------------------

### Delayed and immediate effects as well as chronic effects from short and long-term exposure

Chronic Toxicity	No information available.
------------------	---------------------------

### Numerical measures of toxicity - Product Information

Unknown acute toxicity	100% of the mixture consists of ingredient(s) of unknown toxicity
------------------------	---

## 12. ECOLOGICAL INFORMATION

---

Ecotoxicity	May cause long lasting harmful effects to aquatic life
-------------	--

100% of the mixture consists of component(s) of unknown hazards to the aquatic environment.

Persistence and degradability	No information available.
Bioaccumulation	No information available.
Mobility	No information available.

## 13. DISPOSAL CONSIDERATIONS

---

Disposal of wastes	Disposal should be in accordance with applicable regional, national and local laws and regulations.
Contaminated packaging	Do not reuse container.

---

## 14. TRANSPORT INFORMATION

---

DOT Not regulated

IMDG Not regulated

IATA Not regulated

## 15. REGULATORY INFORMATION

---

### International Inventories

All of the components in the product are on the following Inventory lists

TSCA (United States): Canada (DSL/NDSL) Europe (EINECS/ELINCS/NLP)

Chemical name	TSCA	DSL	NDSL	EINECS	ELINCS	ENCS	IECSC	KECL	PICCS	AICS
Nile Red	X	-	X	X	-	-	-	-	-	-

X - Listed

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

### US Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

#### **SARA 311/312 Hazard Categories**

Acute health hazard No  
Chronic Health Hazard No  
Fire hazard No  
Sudden release of pressure hazard No  
Reactive hazard No

#### **CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

### US State Regulations

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

#### **U.S. State Right-to-Know Regulations**

This product does not contain any substances regulated by state right-to-know regulations

---

## 16. OTHER INFORMATION

---

Revision note No information available

### **Disclaimer**

---

The information provided in this Material Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.