



# SAFETY DATA SHEET

Revision date 01-Jan-2021

Version 1.0

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

### Product identifier

Product Name Tropaeolin OO  
Product Code S10139

### Recommended use of the chemical and restrictions on use

For research use only. Not intended for diagnostic or therapeutic use.

### Details of the supplier of the safety data sheet

DAWN SCIENTIFIC  
121 Liberty Street, Metuchen, NJ  
1-800-DAWN-SCI | 732-902-6300

### Emergency telephone number

Chemtrec  
1.800.424.9300 (Within USA)  
+1.703.527.3887 (Outside USA)

## 2. HAZARDS IDENTIFICATION

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

### Classification

Skin corrosion/irritation Category 2  
Serious eye damage/eye irritation Category 2A  
Specific target organ toxicity (single exposure) Category 3

### Label elements

Signal word Warning  
Hazard statements Causes skin irritation  
Causes serious eye irritation  
May cause respiratory irritation

Symbols/Pictograms



Precautionary Statements - Prevention

Wash face, hands and any exposed skin thoroughly after handling  
Wear protective gloves/protective clothing/eye protection/face protection  
Avoid breathing dust/fume/gas/mist/vapors/spray

Precautionary Statements - Response

Use only outdoors or in a well-ventilated area  
IF exposed or concerned: Get medical advice/attention  
IF IN EYES: Rinse cautiously with water for several minutes.  
Remove contact lenses, if present and easy to do. Continue rinsing  
If eye irritation persists: Get medical advice/attention  
IF ON SKIN: Wash with plenty of soap and water  
If skin irritation occurs: Get medical advice/attention  
Take off contaminated clothing and wash before reuse  
IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing  
Call a POISON CENTER or doctor/physician if you feel unwell

Precautionary Statements - Storage

Store in a well-ventilated place. Keep container tightly closed  
Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

### Hazards not otherwise classified (HNOC)

Hazards not otherwise classified (HNOC)

Not applicable

### Other Information

Unknown acute toxicity

100% of the mixture consists of ingredient(s) of unknown toxicity.

### NFPA

Health hazards	2
Flammability	0
Stability	0
Physical and chemical properties	-



### HMIS

Health hazards	2
Flammability	0
Physical hazards	0
Personal protection	-

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

CAS No	554-73-4
Molecular Weight	375.38
Formula	C <sub>18</sub> H <sub>14</sub> N <sub>3</sub> NaO <sub>3</sub> S

Chemical Name	CAS No	Weight %	Oral LD50	Dermal LD50	Inhalation LC50
Tropaeolin OO	554-73-4	>98	-	-	-

## 4. FIRST AID MEASURES

### First Aid Measures

General advice	Consult a physician if necessary. Remove to fresh air.
Eye contact	Wash with plenty of water.
Skin Contact	Wash skin with soap and water.
Inhalation	Remove to fresh air. If breathing is difficult, give oxygen. If not breathing, give artificial respiration.
Ingestion	Never give anything by mouth to an unconscious person. Clean mouth with water.

### Most important symptoms and effects, both acute and delayed

Symptoms No information available.

### Indication of any immediate medical attention and special treatment needed

Note to physicians Treat symptomatically.

## 5. FIRE-FIGHTING MEASURES

### Suitable Extinguishing Media

Suitable Extinguishing Media Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media None.

### Specific hazards arising from the chemical

Specific hazards arising from the chemical No information available.

Hazardous combustion products Hydrogen sulfide. Carbon oxides. Nitrogen oxides (NO<sub>x</sub>).

---

### Explosion data

Sensitivity to Mechanical Impact No information available.  
Sensitivity to Static Discharge No information available.

### Protective equipment and precautions for firefighters

Protective equipment and precautions for firefighters As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

---

## 6. ACCIDENTAL RELEASE MEASURES

### Personal precautions, protective equipment and emergency procedures

Personal precautions Ensure adequate ventilation, especially in confined areas.

### Environmental precautions

Environmental precautions See Section 12 for additional Ecological Information.

### Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so.  
Methods for cleaning up Use personal protective equipment as required. Cover powder spill with plastic sheet or tarp to minimize spreading and keep powder dry. Take up mechanically, placing in appropriate containers for disposal. Avoid creating dust. Clean contaminated surface thoroughly.

---

## 7. HANDLING AND STORAGE

### Precautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice. Corrosive hazard. Wear protective gloves/clothing and eye/face protection.

### Conditions for safe storage, including any incompatibilities

Storage Conditions Keep containers tightly closed in a dry, cool and well-ventilated place. Store at room temperature.  
Incompatible materials None known based on information supplied.

---

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### Control parameters

Exposure Guidelines This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

### Appropriate engineering controls

Engineering Controls Showers  
Eyewash stations  
Ventilation systems

### Individual protection measures, such as personal protective equipment

Eye/face protection Wear safety glasses with side shields (or goggles).  
Skin and Body Protection Wear protective gloves and protective clothing.  
Respiratory protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.  
General Hygiene Considerations Handle in accordance with good industrial hygiene and safety practice.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

---

Physical State	Solid
Appearance	No information available
Odor	No information available

### **Property**

pH	No information available
Melting point/freezing point	No information available
Boiling point	No information available
Flash point	No information available
Density	No information available
Evaporation rate	No information available
Upper flammability limits	No information available
Lower flammability limit	No information available
Vapor pressure	No information available
Vapor density	No information available
Specific gravity	No information available
Water solubility	No information available
Solubility in other solvents	No information available
Partition coefficient	No information available
Autoignition temperature	No information available
Decomposition temperature	No information available
Kinematic viscosity	No information available
Explosive properties	No information available
Oxidizing properties	No information available

### **Values**

---

## 10. STABILITY AND REACTIVITY

---

Reactivity	Not applicable
Chemical stability	Stable under recommended storage conditions.
Possibility of Hazardous Reactions	None under normal processing.
Hazardous polymerization	No information available.
Conditions to avoid	Extremes of temperature and direct sunlight.
Incompatible materials	Strong oxidizing agents.
Hazardous Decomposition Products	Hydrogen sulfide. Carbon oxides. Nitrogen oxides (NOx).

---

## 11. TOXICOLOGICAL INFORMATION

---

### **Information on likely routes of exposure**

Inhalation	No data available.
Eye contact	No data available.
Skin Contact	No data available.
Ingestion	No data available.

### **Information on toxicological effects**

Symptoms	No information available.
----------	---------------------------

### **Delayed and immediate effects as well as chronic effects from short and long-term exposure**

Chronic Toxicity	No information available.
------------------	---------------------------

### **Numerical measures of toxicity - Product Information**

Unknown acute toxicity	100% of the mixture consists of ingredient(s) of unknown toxicity
------------------------	---

---

## 12. ECOLOGICAL INFORMATION

---

Ecotoxicity	May cause long lasting harmful effects to aquatic life
100% of the mixture consists of components(s) of unknown hazards to the aquatic environment.	
Persistence and degradability	No information available.
Bioaccumulation	No information available.
Mobility	No information available.

## 13. DISPOSAL CONSIDERATIONS

---

Disposal of wastes	Disposal should be in accordance with applicable regional, national and local laws and regulations.
Contaminated packaging	Do not reuse container.

## 14. TRANSPORT INFORMATION

---

DOT	Not regulated
IMDG	Not regulated
IATA	Not regulated

## 15. REGULATORY INFORMATION

---

### International Inventories

**All of the components in the product are on the following Inventory lists**

TSCA (United States): Canada (DSL/NDSL) Europe (EINECS/ELINCS/NLP) Australia (AICS) South Korea (KECL): China (IECSC)

Chemical Name	TSCA	DSL	NDSL	EINECS	ELINCS	ENCS	IECSC	KECL	PICCS	AICS
Tropaeolin OO	X	-	X	X	-	-	X	X	-	X

X - Listed

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

### US Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

#### **SARA 311/312 Hazard Categories**

Acute health hazard No

