



# SAFETY DATA SHEET

Revision date 01-Jan-2021  
Version 1.0

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

### Product identifier

Product Name Pararosaniline Base  
Product Code S10142

### Recommended use of the chemical and restrictions on use

For research use only. Not intended for diagnostic or therapeutic use.

### Details of the supplier of the safety data sheet

DAWN SCIENTIFIC  
121 Liberty Street, Metuchen, NJ  
1-800-DAWN-SCI | 732-902-6300

### Emergency telephone number

Chemtrec  
1.800.424.9300 (Within USA)  
+1.703.527.3887 (Outside USA)

## 2. HAZARDS IDENTIFICATION

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

### Classification

Carcinogenicity Category 2

### Label elements

Signal word Warning  
Hazard statements Suspected of causing cancer  
Symbols/Pictograms



Precautionary Statements - Prevention

Obtain special instructions before use  
Do not handle until all safety precautions have been read and understood

Precautionary Statements - Response

Use personal protective equipment as required  
IF exposed or concerned: Get medical advice/attention

Precautionary Statements - Storage

Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

### Hazards not otherwise classified (HNOC)


Hazards not otherwise classified (HNOC) Not applicable

### Other Information

Unknown acute toxicity 100% of the mixture consists of ingredient(s) of unknown toxicity.

### NFPA

Health hazards -  
Flammability -  
Stability -  
Physical and chemical properties -



### HMIS

Health hazards -  
Flammability -  
Physical hazards -  
Personal protection -  
\* = Chronic Health Hazard

---

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

---

CAS No 467-62-9  
Molecular Weight 305.37  
Formula C<sub>19</sub>H<sub>19</sub>N<sub>3</sub>O

Chemical name	CAS No.	Weight-%	Oral LD50	Dermal LD50	Inhalation LC50
Pararosaniline Base	467-62-9	>98	-	-	-

---

### 4. FIRST AID MEASURES

---

#### First Aid Measures

General advice Consult a physician if necessary. Remove to fresh air.  
Eye contact Wash with plenty of water.  
Skin Contact Wash skin with soap and water.  
Inhalation Remove to fresh air. If breathing is difficult, give oxygen. If not breathing, give artificial respiration.  
Ingestion Never give anything by mouth to an unconscious person. Clean mouth with water.

#### Most important symptoms and effects, both acute and delayed

Symptoms No information available.

#### Indication of any immediate medical attention and special treatment needed

Note to physicians Treat symptomatically.

---

### 5. FIRE-FIGHTING MEASURES

---

#### Suitable Extinguishing Media

Suitable Extinguishing Media Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

#### Specific hazards arising from the chemical

Specific hazards arising from the chemical No information available.

Hazardous combustion products Carbon oxides. Nitrogen oxides (NO<sub>x</sub>).

#### Explosion data

Sensitivity to Mechanical Impact No information available.

Sensitivity to Static Discharge No information available.

#### Protective equipment and precautions for firefighters

Protective equipment and precautions for firefighters As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

---

### 6. ACCIDENTAL RELEASE MEASURES

---

#### Personal precautions, protective equipment and emergency procedures

Personal precautions Ensure adequate ventilation, especially in confined areas.

#### Environmental precautions

Environmental precautions See Section 12 for additional Ecological Information.

---

### Methods and material for containment and cleaning up

Methods for containment	Prevent further leakage or spillage if safe to do so.
Methods for cleaning up	Use personal protective equipment as required. Cover powder spill with plastic sheet or tarp to minimize spreading and keep powder dry. Take up mechanically, placing in appropriate containers for disposal. Avoid creating dust. Clean contaminated surface thoroughly.

## 7. HANDLING AND STORAGE

---

### Precautions for safe handling

Advice on safe handling	Handle in accordance with good industrial hygiene and safety practice.
-------------------------	--

### Conditions for safe storage, including any incompatibilities

Storage Conditions	Keep containers tightly closed in a dry, cool and well-ventilated place. Store at room temperature.
Incompatible materials	None known based on information supplied.

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

---

### Control parameters

Exposure Guidelines	This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.
---------------------	---

### Appropriate engineering controls

Engineering Controls	Showers Eyewash stations Ventilation systems
----------------------	--

### Individual protection measures, such as personal protective equipment

Eye/face protection	Wear safety glasses with side shields (or goggles).
Skin and Body Protection	Wear protective gloves and protective clothing.
Respiratory protection	If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.
General Hygiene Considerations	Handle in accordance with good industrial hygiene and safety practice.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

---

Physical State	Solid
Appearance	powder
Odor	Odorless

### Property

<u>Values</u>	
pH	No information available
Melting point/freezing point	205 °C
Boiling point	No information available
Flash point	No information available
Liquid Density	No information available
Evaporation rate	No information available
Upper flammability limits	No information available
Lower flammability limit	No information available
Vapor pressure	No information available

---

Vapor density	1
Specific gravity	1.33
Water solubility	No information available
Solubility in other solvents	No information available
Partition coefficient	1.63
Autoignition temperature	No information available
Decomposition temperature	205 ° C
Kinematic viscosity	No information available
Explosive properties	No information available
Oxidizing properties	No information available

## 10. STABILITY AND REACTIVITY

---

Reactivity	Not applicable
Chemical stability	Stable under recommended storage conditions.
Possibility of Hazardous Reactions	None under normal processing.
Hazardous polymerization	No information available.
Conditions to avoid	Extremes of temperature and direct sunlight.
Incompatible materials	Strong oxidizing agents.
Hazardous Decomposition Products	Carbon oxides. Nitrogen oxides (NOx).

## 11. TOXICOLOGICAL INFORMATION

---

### Information on likely routes of exposure

Inhalation	No data available.
Eye contact	No data available.
Skin Contact	No data available.
Ingestion	No data available.

### Information on toxicological effects

Symptoms	No information available.
----------	---------------------------

### Delayed and immediate effects as well as chronic effects from short and long-term exposure

Chronic Toxicity	No information available.
------------------	---------------------------

### Numerical measures of toxicity - Product Information

Unknown acute toxicity	100% of the mixture consists of ingredient(s) of unknown toxicity
------------------------	---

## 12. ECOLOGICAL INFORMATION

---

Ecotoxicity	May cause long lasting harmful effects to aquatic life
100% of the mixture consists of components(s) of unknown hazards to the aquatic environment.	
Persistence and degradability	No information available.
Bioaccumulation	No information available.
Mobility	No information available.

## 13. DISPOSAL CONSIDERATIONS

---

Disposal of wastes Disposal should be in accordance with applicable regional, national and local laws and regulations.  
Contaminated packaging Do not reuse container.

## 14. TRANSPORT INFORMATION

---

**DOT** Not regulated  
**IMDG** Not regulated  
**IATA** Not regulated

## 15. REGULATORY INFORMATION

---

### **International Inventories**

**All of the components in the product are on the following Inventory lists**

TSCA (United States): Canada (DSL/NDSL) Europe (EINECS/ELINCS/NLP) South Korea (KECL): Philippines (PICCS)

Chemical name	TSCA	DSL	NDSL	EINECS	ELINCS	ENCS	IECSC	KECL	PICCS	AICS
Pararosaniline Base	X	-	X	X	-	-	-	X	X	-

X - Listed

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

### **US Federal Regulations**

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

#### **SARA 311/312 Hazard Categories**

Acute health hazard No  
Chronic Health Hazard Yes  
Fire hazard No  
Sudden release of pressure hazard No  
Reactive hazard No

#### **CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

### **US State Regulations**

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

#### **U.S. State Right-to-Know Regulations**

This product does not contain any substances regulated by state right-to-know regulations

---

## 16. OTHER INFORMATION

---

Revision note

No information available

### **Disclaimer**

The information provided in this Material Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.