



Tri Stains

# SAFETY DATA SHEET

Revision date 01-Jan-2021  
Version 1.0

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

### Product identifier

Product Name Patent Blue VF  
Product Code S10147

### Recommended use of the chemical and restrictions on use

For research use only. Not intended for diagnostic or therapeutic use.

### Details of the supplier of the safety data sheet

DAWN SCIENTIFIC  
121 Liberty Street, Metuchen, NJ  
1-800-DAWN-SCI | 732-902-6300

### Emergency telephone number

Chemtrec  
1.800.424.9300 (Within USA)  
+1.703.527.3887 (Outside USA)

## 2. HAZARDS IDENTIFICATION

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

### Classification

Serious eye damage/eye irritation  
Specific target organ toxicity (single exposure)

Category 2A  
Category 3

### Label elements

Signal word  
Hazard statements

Warning  
Causes serious eye irritation  
May cause respiratory irritation

Symbols/Pictograms



Precautionary Statements - Prevention

Wash face, hands and any exposed skin thoroughly after handling  
Wear protective gloves/protective clothing/eye protection/face protection  
Avoid breathing dust/fume/gas/mist/vapors/spray

Precautionary Statements - Response

Use only outdoors or in a well-ventilated area  
IF exposed or concerned: Get medical advice/attention  
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing  
If eye irritation persists: Get medical advice/attention  
IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

Precautionary Statements - Storage

Call a POISON CENTER or doctor/physician if you feel unwell  
Store in a well-ventilated place. Keep container tightly closed  
Store locked up

Precautionary Statements - Disposal


Dispose of contents/container to an approved waste disposal plant

### Hazards not otherwise classified (HNOC)

Hazards not otherwise classified (HNOC)

Not applicable

### Other Information

<b>NFPA</b>	Health hazards	0		<b>HMIS</b>	Health hazards	0	
	Flammability	0				Flammability	0
	Stability	0				Physical hazards	0
	Physical and chemical properties	-				Personal protection	-

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

CAS No 129-17-9  
Molecular Weight 566.67  
Formula  $C_{27}H_{31}N_2NaO_6S_2$

Chemical Name	CAS No	Weight %	Oral LD50	Dermal LD50	Inhalation LC50
Patent Blue VF	129-17-9	>98	-	-	-

## 4. FIRST AID MEASURES

### First Aid Measures

General advice	If symptoms persist, call a physician.
Eye contact	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. If symptoms persist, call a physician. Immediately flush with plenty of water. After initial flushing, remove any contact lenses and continue flushing for at least 15 minutes. Keep eye wide open while rinsing.
Skin Contact	Wash off immediately with plenty of water for at least 15 minutes. If skin irritation persists, call a physician. Wash contaminated clothing before reuse. Wash off immediately with plenty of water. Immediate medical attention is not required. Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes.
Inhalation	Remove to fresh air. If symptoms persist, call a physician. Artificial respiration and/or oxygen may be necessary. Call a physician. Immediate medical attention is not required. Move to fresh air in case of accidental inhalation of vapors.
Ingestion	Immediate medical attention is not required. Rinse mouth. Drink plenty of water. Do NOT induce vomiting. Clean mouth with water and drink afterwards plenty of water. Never give anything by mouth to an unconscious person. Call a physician.
Self-protection of the first aider	Use personal protective equipment as required.

### Most important symptoms and effects, both acute and delayed

Symptoms No information available.

### Indication of any immediate medical attention and special treatment needed

Note to physicians Treat symptomatically.

## 5. FIRE-FIGHTING MEASURES

### Suitable Extinguishing Media

Suitable Extinguishing Media Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media None.

### Specific hazards arising from the chemical

---

Specific hazards arising from the chemical	No information available.
Hazardous combustion products	Hydrogen sulfide. Carbon oxides. Nitrogen oxides (NOx).

#### **Explosion data**

Sensitivity to Mechanical Impact	No information available.
Sensitivity to Static Discharge	No information available.

#### **Protective equipment and precautions for firefighters**

Protective equipment and precautions for firefighters	As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.
---	--

---

## **6. ACCIDENTAL RELEASE MEASURES**

#### **Personal precautions, protective equipment and emergency procedures**

Personal precautions	Use personal protective equipment as required. Avoid contact with eyes and skin. Evacuate personnel to safe areas. Keep people away from and upwind of spill/leak.
----------------------	--

#### **Environmental precautions**

Environmental precautions	Prevent entry into waterways, sewers, basements or confined areas. Do not flush into surface water or sanitary sewer system. Prevent further leakage or spillage if safe to do so. Prevent product from entering drains. See Section 12 for additional Ecological Information.
---------------------------	--

#### **Methods and material for containment and cleaning up**

Methods for containment	Prevent further leakage or spillage if safe to do so.
Methods for cleaning up	Use personal protective equipment as required. Cover powder spill with plastic sheet or tarp to minimize spreading and keep powder dry. Take up mechanically, placing in appropriate containers for disposal. Avoid creating dust. Clean contaminated surface thoroughly. Soak up with inert absorbent material. Dam up. Pick up and transfer to properly labeled containers. Sweep up and shovel into suitable containers for disposal. After cleaning, flush away traces with water. Take precautionary measures against static discharges.

---

## **7. HANDLING AND STORAGE**

#### **Precautions for safe handling**

Advice on safe handling	Ensure adequate ventilation, especially in confined areas. Use personal protective equipment as required. Use with local exhaust ventilation. Do not breathe dust/fume/gas/mist/vapors/spray. Noxious vapor/odor.
-------------------------	---

#### **Conditions for safe storage, including any incompatibilities**

Storage Conditions	Keep container tightly closed. Keep out of the reach of children. Keep containers tightly closed in a cool, well-ventilated place. Keep in properly labeled containers. Store at room temperature.
Incompatible materials	None known based on information supplied.

---

## **8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

#### **Control parameters**

Exposure Guidelines	This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.
---------------------	---

#### **Appropriate engineering controls**

---

Engineering Controls

Ensure adequate ventilation, especially in confined areas

**Individual protection measures, such as personal protective equipment**

Eye/face protection

Tight sealing safety goggles. Face protection shield.

Skin and Body Protection

Wear protective gloves and protective clothing.

Respiratory protection

If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

General Hygiene Considerations

When using do not eat, drink or smoke. Wash contaminated clothing before reuse. Regular cleaning of equipment, work area and clothing is recommended.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

Physical State

Solid

Appearance

powder

Odor

May be unpleasant

**Property**

**Values**

pH

No information available

Melting point/freezing point

300 °C

Boiling point

No information available

Flash point

No information available

Density

No information available

Evaporation rate

No information available

Upper flammability limits

No information available

Lower flammability limit

No information available

Vapor pressure

No information available

Vapor density

No information available

Specific gravity

No information available

Water solubility

No information available

Solubility in other solvents

No information available

Partition coefficient

0.57

Autoignition temperature

No information available

Decomposition temperature

No information available

Kinematic viscosity

No information available

Explosive properties

No information available

Oxidizing properties

No information available

---

## 10. STABILITY AND REACTIVITY

Reactivity

Not applicable

Chemical stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Hazardous polymerization

No information available.

Conditions to avoid

Extremes of temperature and direct sunlight.

Incompatible materials

Strong oxidizing agents.

Hazardous Decomposition Products

Hydrogen sulfide. Carbon oxides. Nitrogen oxides (NOx).

---

## 11. TOXICOLOGICAL INFORMATION

**Information on likely routes of exposure**

Inhalation

No data available.

Eye contact

No data available.

Skin Contact

No data available.

---

Ingestion No data available.

**Information on toxicological effects**

Symptoms No information available.

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

Chronic Toxicity Avoid repeated exposure.

Carcinogenicity The table below indicates whether each agency has listed any ingredient as a carcinogen.

Chemical Name	ACGIH	IARC	NTP	OSHA
Patent Blue VF 129-17-9	-	Group 3	-	-

*IARC (International Agency for Research on Cancer) Not classifiable as a human carcinogen*

**Numerical measures of toxicity - Product Information**

Unknown acute toxicity No information available

---

**12. ECOLOGICAL INFORMATION**

Ecotoxicity May cause long lasting harmful effects to aquatic life

100% of the mixture consists of component(s) of unknown hazards to the aquatic environment.

Persistence and degradability No information available.

Bioaccumulation No information available.

Mobility No information available.

---

**13. DISPOSAL CONSIDERATIONS**

Disposal of wastes Disposal should be in accordance with applicable regional, national and local laws and regulations.

Contaminated packaging Do not reuse container.

Other Information Waste codes should be assigned by the user based on the application for which the product was used.

---

**14. TRANSPORT INFORMATION**

DOT Not regulated

IMDG Not regulated

IATA Not regulated

---

**15. REGULATORY INFORMATION**

**International Inventories**

**All of the components in the product are on the following Inventory lists**

TSCA (United States): Canada (DSL/NDL) Europe (EINECS/ELINCS/NLP) Australia (AICS) South Korea (KECL): China (IECSC)  
ENCS (Japan): Philippines (PICCS)

Chemical Name	TSCA	DSL	NDSL	EINECS	ELINCS	ENCS	IECSC	KECL	PICCS	AICS
Patent Blue VF	X	X	-	X	-	X	X	X	X	X

X - Listed

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

## US Federal Regulations

### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

### **SARA 311/312 Hazard Categories**

Acute health hazard	Yes
Chronic Health Hazard	No
Fire hazard	No
Sudden release of pressure hazard	No
Reactive hazard	No

### **CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

## US State Regulations

### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

### **U.S. State Right-to-Know Regulations**

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Patent Blue VF 129-17-9	X	Not Listed	Not Listed

## **16. OTHER INFORMATION**

Revision note

No information available

### **Disclaimer**

The information provided in this Material Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.