

Preparation Date: 01/01/2019

Revision Date: N/A

Revision Number: N/A

1. IDENTIFICATION**Product identifier****Product code:** C7470**Product Name:** SODIUM CITRATE, DIHYDRATE, CRYSTAL, REAGENT, ACS**Other means of identification****Synonyms:** Trisodium citrate dihydrate
1,2,3-Propanetricarboxylic acid, 2-hydroxy-, trisodium salt, dihydrate**CAS #:** 6132-04-3**RTECS #** TZ8595800 (Dihydrate)
GE8300000 (Anhydrous)**CI#:** Not available**Recommended use of the chemical and restrictions on use****Recommended use:** No information available.**Uses advised against** No information available**Supplier:****Dawn Scientific Inc**

121 Liberty Street, Metuchen, NJ, 08840

Tel : 732-902-6300 | Fax : 973-802-1005

sales@dawnscientific.com | www.dawnscientific.com

Emergency telephone number Chemtrec 1-800-424-9300**2. HAZARDS IDENTIFICATION****Classification**

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS)

Label elements**Not classified****Hazards not otherwise classified (HNOC)**

Not Applicable

Other hazards

Not available

3. COMPOSITION/INFORMATION ON INGREDIENTS

Component	CAS No	Weight-%
Sodium Citrate, Dihydrate	6132-04-3	100

4. FIRST AID MEASURES

First aid measures

General Advice:	National Capital Poison Center in the United States can provide assistance if you have a poison emergency and need to talk to a poison specialist. Call 1-800-222-1222.
Skin Contact:	Wash off immediately with soap and plenty of water removing all contaminated clothing and shoes. Get medical attention if irritation develops. Consult a physician if necessary.
Eye Contact:	Flush eyes with water for 15 minutes. Get medical attention if irritation occurs. If symptoms persist, call a physician.
Inhalation:	Move to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention.
Ingestion:	Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person. Consult a physician if necessary.

Most important symptoms and effects, both acute and delayed

Symptoms	May cause eye/skin irritation May cause digestive (gastrointestinal) tract irritation May cause gastrointestinal disturbances
-----------------	---

Indication of any immediate medical attention and special treatment needed

Notes to Physician:	Treat symptomatically.
----------------------------	------------------------

Protection of first-aiders

First-Aid Providers: Avoid exposure to blood or body fluids. Wear gloves and other necessary protective clothing. Dispose of contaminated clothing and equipment as bio-hazardous waste.

5. FIRE-FIGHTING MEASURES

Extinguishing Media

Suitable Extinguishing Media:	Carbon dioxide (CO2). Dry chemical. Water spray mist or foam.
--------------------------------------	---

Unsuitable Extinguishing Media:	No information available.
--	---------------------------

Specific hazards arising from the chemical

Hazardous combustion products	Carbon Monoxide, Carbon Dioxide. Sodium oxides.
Specific hazards	May be combustible at high temperatures.

Special Protective Actions for Firefighters

Specific Methods:

No information available

Special Protective Equipment for Firefighters:

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal Precautions:

Ensure adequate ventilation. Use personal protective equipment. Avoid contact with skin, eyes and clothing. Avoid dust formation. Remove all sources of ignition.

Environmental precautions

Prevent further leakage or spillage if safe to do so. Prevent entry into waterways, sewers, basements or confined areas.

Methods and material for containment and cleaning up

Methods for containment

Stop leak if you can do it without risk. Cover with plastic sheet to prevent spreading.

Methods for cleaning up

Sweep up and shovel into suitable containers for disposal. Clean contaminated surface thoroughly.

7. HANDLING AND STORAGE

Precautions for safe handling

Technical Measures/Precautions:

Provide sufficient air exchange and/or exhaust in work rooms. Avoid dust formation. Keep away from incompatible materials.

Safe Handling Advice:

Wear personal protective equipment. Avoid contact with skin, eyes and clothing. Keep away from heat and sources of ignition. Avoid dust formation. Do not ingest. Do not breathe dust. Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Technical Measures/Storage Conditions:

Keep container tightly closed in a dry and well-ventilated place. Store at room temperature in the original container. Store away from incompatible materials.

Incompatible Materials:

Strong oxidizing agents

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

National occupational exposure limits

United States

Component	CAS No	OSHA	NIOSH	ACGIH	AIHA WEEL
Sodium Citrate,	6132-04-3	None	None	None	None

Dihydrate					
-----------	--	--	--	--	--

Canada

Component	CAS No	Canada - Alberta	Canada - British Columbia	Canada - Ontario	Canada - Quebec
Sodium Citrate, Dihydrate	6132-04-3	None	None	None	None

Australia and Mexico

Component	CAS No	Australia	Mexico
Sodium Citrate, Dihydrate	6132-04-3	None	None

Appropriate engineering controls

Engineering measures to reduce exposure:

Ensure adequate ventilation. Use process enclosures, local exhaust ventilation, or other engineering controls to keep airborne levels below recommended exposure limits. If user operations generate dust, fume or mist, use ventilation to keep exposure to airborne contaminants below the exposure limit.

Individual protection measures, such as personal protective equipment

Personal Protective Equipment

Eye protection:	Safety glasses with side-shields.
Skin and body protection:	Long sleeved clothing Chemical resistant apron Gloves
Respiratory protection:	Effective dust mask. Use a dust respirator under conditions where exposure to the substance is apparent (e.g. generation of high concentration of dust (dust clouds) , inadequate ventilation, development of respiratory tract irritation), and engineering controls are not feasible. Be sure to use an approved/certified respirator or equivalent.
Hygiene measures:	Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product When using, do not eat, drink or smoke.

9. PHYSICAL AND CHEMICAL PROPERTIES

Physical state: Solid	Appearance: Granular.	Color: White.
Odor: Odorless.	Taste Cool. Saline.	Formula C ₆ H ₅ -Na-3O ₇ .2H ₂ O
Molecular/Formula weight (g/mole): 294.1	Flammability (solid, gas) no data available	Flashpoint (°C/°F): No information available
Flash Point Tested according to: Not available	Autoignition Temperature (°C/°F): No information available	Lower Explosion Limit (%): No information available
Upper Explosion Limit (%):		

No information available

Melting point/range(°C/°F):
313.9 °C/572 °F

Decomposition temperature(°C/°F):
Decomposes at red heat Becomes anhydrous at 150 °C

Boiling point/range(°C/°F):
No information available

Bulk density:
No information available

Density (g/cm3):
No information available

Specific gravity:
1.665

pH
No information available

Vapor pressure @ 20°C (kPa):
No information available

Evaporation rate:
No information available

Vapor density:
No information available

VOC content (g/L):
No information available

Odor threshold (ppm):
No information available

Partition coefficient (n-octanol/water):
No information available

Viscosity:
No information available

Miscibility:
No information available

Solubility:
Freely soluble in water
Insoluble in Alcohol

10. STABILITY AND REACTIVITY

Reactivity

Reactive with oxidizing agents

Chemical stability

Stability: Stable under recommended storage conditions.

Possibility of Hazardous Reactions: Hazardous polymerization does not occur

Conditions to avoid: Heat. Avoid dust formation. Incompatible materials.

Incompatible Materials: Strong oxidizing agents

Hazardous decomposition products: Carbon monoxide. Carbon dioxide. Sodium oxides.

Other Information

Corrosivity: No information available

Special Remarks on Corrosivity: No information available

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Principal Routes of Exposure:
Ingestion. Inhalation.

Acute Toxicity

Component Information

Sodium Citrate, Dihydrate

CAS No 6132-04-3

LD50/oral/rat = No information available
LD50/oral/mouse = 5400 mg/kg (ECHA - Registration Dossier)
LD50/dermal/rabbit = No information available
LD50/dermal/rat = >2000 mg/kg (ECHA - Registration Dossier)
LC50/inhalation/rat = No information available
LC50/inhalation/mouse = No information available
Other LD50 or LC50 information = No information available

Product Information

LD50/oral/rat =
Value - Acute Toxicity = No information available

LD50/oral/mouse =
Value - Acute Tox = 5400 mg/kg

LD50/dermal/rabbit
Value - Acute Toxicity = No information available

LD50/dermal/rat
VALUE - Acute Tox = >2000 mg/kg

LC50/inhalation/rat
VALUE-Vapor = No information available
VALUE-Gas = No information available
VALUE-Dust/Mist = No information available

LC50/Inhalation/mouse
VALUE-Vapor = No information available
VALUE - Gas = No information available
VALUE - Dust/Mist = No information available

Symptoms

Skin Contact:	May cause skin irritation.
Eye Contact:	May cause eye irritation.
Inhalation	May cause irritation of respiratory tract.
Ingestion	Low hazard. May cause digestive (gastrointestinal) tract irritation. May cause gastrointestinal disturbances. May cause alkalosis. It may depress the heart by decreasing the ionized calcium level of the blood. It may affect behavior/central nervous system (tetany, convulsions).
Aspiration hazard	No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Chronic Toxicity	No information available.
Sensitization:	No information available.
Mutagenic Effects:	No information available

Carcinogenic effects: Not considered carcinogenic.

Component	CAS No	IARC	ACGIH - Carcinogens	NTP	OSHA HCS - Carcinogens	Australia - Notifiable Carcinogenic Substances	Australia - Prohibited Carcinogenic Substances
Sodium Citrate, Dihydrate	6132-04-3	Not listed	Not listed	Not listed	Not listed	Not listed	Not listed

ACGIH (American Conference of Governmental Industrial Hygienists)

IARC (International Agency for Research on Cancer)

NTP (National Toxicology Program)

OSHA (Occupational Safety and Health Administration of the US Department of Labor)

Reproductive toxicity No data is available

Reproductive Effects: No information available

Developmental Effects: No information available

Teratogenic Effects: No information available

Specific Target Organ Toxicity

STOT - single exposure No information available.

STOT - repeated exposure No information available.

Target Organs: No information available.

12. ECOLOGICAL INFORMATION

Ecotoxicity

Ecotoxicity effects: No data available.

Persistence and degradability: No information available

Bioaccumulative potential: No information available.

Mobility in soil No information available

Other adverse effects No information available.

13. DISPOSAL CONSIDERATIONS

Disposal Methods

Waste from residues / unused products:

Waste must be disposed of in accordance with Federal, State and Local regulation.

Contaminated packaging:

Empty containers should be taken for local recycling, recovery or waste disposal

Component	CAS No	RCRA - F Series Wastes	RCRA - K Series Wastes	RCRA - P Series Wastes	RCRA - U Series Wastes
Sodium Citrate, Dihydrate	6132-04-3	None	None	None	None

14. TRANSPORT INFORMATION

DOT

UN-No:	Not Regulated
Proper Shipping Name:	No information available
Hazard Class	No information available
Subsidiary Class	No information available
Packing group:	No information available
Emergency Response Guide Number	No information available
Marine Pollutant	No data available
DOT RQ (lbs):	No information available
Special Provisions	No Information available
Symbol(s):	No information available
Description:	No information available

TDG (Canada)

UN-No:	Not Regulated
Proper Shipping Name:	No information available
Hazard Class	No information available
Subsidiary Risk:	No information available
Packing Group:	No information available
Marine Pollutant	No Information available
Description:	No information available

ADR

UN Number	Not regulated
Proper Shipping Name:	No information available
Transport hazard class(es)	No information available
Packing group	No information available
Subsidiary Risk:	No information available

IMDG

UN-No:	Not Regulated
Proper Shipping Name:	No information available
Hazard Class:	No information available
Subsidiary Risk:	No information available
Packing Group:	No information available
Marine Pollutant	No information available

RID

UN Number	Not Regulated
Proper Shipping Name:	No information available
Transport hazard class(es)	No information available
Subsidiary Risk:	No information available
Packing group	No information available

ICAO (air)

UN-No:	Not Regulated
Proper Shipping Name:	No information available
Hazard Class	No information available
Subsidiary Risk:	No information available
Packing Group:	No information available

IATA

UN Number	Not Regulated
Proper Shipping Name:	No information available

Transport hazard class(es) No information available
Subsidiary Risk: No information available
Packing group No information available
Precautionary Statements - Response No information available
Special Provisions No information available

15. REGULATORY INFORMATION

International Inventories

Component	CAS No	U.S. TSCA	KOREA KECL	Philippines (PICCS)	Japan ENCS	China IECSC	Australia (AICS)	EINECS-No.
Sodium Citrate, Dihydrate	6132-04-3	Not listed. It is exempt from listing on the TSCA 8(b) Inventory since it is a hydrate. However, the anhydrous form (CAS number 68-04-2) is listed as ACTIVE on the TSCA 8(b) Inventory.	Not present	Present	Not present	Present	Present	Not present

U.S. Regulations

California Prop. 65: Safe Drinking Water and Toxic Enforcement Act of 1986.

Chemicals Known to the State of California to Cause Cancer:

This product does not contain a chemical requiring a warning under California Prop. 65. (See table below)

Chemicals Known to the State of California to Cause Reproductive Toxicity:

This product does not contain a chemical requiring a warning under California Prop. 65. (See table below)

Component	CAS No	Carcinogen	Developmental Toxicity	Male Reproductive Toxicity	Female Reproductive Toxicity:
Sodium Citrate, Dihydrate	6132-04-3	Not Listed	Not Listed	Not Listed	Not Listed

CERCLA/SARA

Component	CAS No	CERCLA - Hazardous Substances and their Reportable Quantities	Section 302 Extremely Hazardous Substances and TPQs	Section 302 Extremely Hazardous Substances and RQs	Section 313 - Chemical Category	Section 313 - Reporting de minimis
Sodium Citrate, Dihydrate	6132-04-3	None	None	None	None	None

U.S. TSCA

Component	CAS No	TSCA Section 5(a)2 - Chemicals With Significant New Use Rules (SNURS)	TSCA 8(d) -Health and Safety Reporting

Sodium Citrate, Dihydrate	6132-04-3	Not Applicable	Not Applicable
---------------------------	-----------	----------------	----------------

Canada

WHMIS 2015 - GHS Classifications

WHMIS 2015 Hazard Classification Information:

Component
Sodium Citrate, Dihydrate
6132-04-3 (100)

WHMIS 2015 Hazard Classification
Combustible Dust - Category 1: May form combustible dust concentrations in air (factors such as combustibility and explosiveness of dusts including composition and shape and size of particles could cause substance to belong to 'Combustible dust' hazard class)

Canada Hazardous Products Regulation This product has been classified according to the hazard criteria of the HPR (Hazardous Products Regulation) and the SDS contains all of the information required by the HPR

DSL/NDSL

Component	CAS No	Canada (DSL)	Canada (NDSL)
Sodium Citrate, Dihydrate	6132-04-3	Not Listed	Not Listed

Component	CAS No	CEPA Schedule I - Toxic Substances
Sodium Citrate, Dihydrate	6132-04-3	Not listed
Component	CAS No	CEPA - 2010 Greenhouse Gases Subject to Mandatory Reporting
Sodium Citrate, Dihydrate	6132-04-3	Not listed

EU Classification

EU GHS - SV - CLP 1272/2008

Component	CAS No	EU GHS - SV - CLP (1272/2008)
Sodium Citrate, Dihydrate	6132-04-3	

EU - CLP (1272/2008)

R-phrase(s)

not determined (not applicable)

S -phrase(s)

none

Component	CAS No	Classification	Concentration Limits:	Safety Phrases
Sodium Citrate, Dihydrate	6132-04-3		No information	

The product is classified in accordance with Annex VI to Directive 67/548/EEC

Indication of danger:

Not dangerous

16. OTHER INFORMATION

Preparation Date:	01/01/2019
Revision date	N/A
Prepared by:	-

Disclaimer:

All chemicals may pose unknown hazards and should be used with caution. This Safety Data Sheet (SDS) applies only to the material as packaged. If this product is combined with other materials, deteriorates, or becomes contaminated, it may pose hazards not mentioned in this SDS. The physical properties reported in this SDS are obtained from the literature and do not constitute product specifications. Information contained herein does not constitute a warranty, whether expressed or implied, as to the safety, merchantability or fitness of the goods for a particular purpose. Dawn Scientific Inc Chemicals & Laboratory Products, assumes no responsibility for results obtained or for incidental or consequential damages, including lost profits, arising from the use of these data. No warranty against infringement of any patent, copyright or trademark is made or implied. It shall be the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use. While this SDS is based on technical data judged to be reliable, Dawn Scientific Inc assumes no responsibility for the completeness or accuracy of the information contained herein.

End of Safety Data Sheet