



# SAFETY DATA SHEET

Revision date 01-Jan-2021  
Version 1.0

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

### Product identifier

Product Name Flavianic acid  
Product Code S10086

### Recommended use of the chemical and restrictions on use

For research use only. Not intended for diagnostic or therapeutic use.

### Details of the supplier of the safety data sheet

DAWN SCIENTIFIC  
121 Liberty Street, Metuchen, NJ  
1-800-DAWN-SCI | 732-902-6300

### Emergency telephone number

Chemtrec  
+1.800.424.9300 (Within USA)  
+1.703.527.3887 (Outside USA)

## 2. HAZARDS IDENTIFICATION

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

### Classification

Skin corrosion/irritation  
Serious eye damage/eye irritation  
Carcinogenicity

Category 1 Sub-category B  
Category 1  
Category 2

### Label elements

Signal word  
Hazard statements

Danger  
Causes severe skin burns and eye damage  
Suspected of causing cancer

Symbols/Pictograms



Precautionary Statements - Prevention

Obtain special instructions before use  
Do not handle until all safety precautions have been read and understood  
Use personal protective equipment as required  
Do not breathe dust/fume/gas/mist/vapors/spray  
Wash face, hands and any exposed skin thoroughly after handling

Precautionary Statements - Response

Immediately call a POISON CENTER or doctor/physician  
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing  
Immediately call a POISON CENTER or doctor/physician  
IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower  
Wash contaminated clothing before reuse  
IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing  
Immediately call a POISON CENTER or doctor/physician  
IF SWALLOWED: Rinse mouth. DO NOT induce vomiting  
Store locked up

Precautionary Statements - Storage

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

### Hazards not otherwise classified (HNOC)

Hazards not otherwise classified (HNOC)

Not applicable

### Other Information

<b>NFPA</b>	Health hazards	3		<b>HMIS</b>	Health hazards	3
	Flammability	0			Flammability	0
	Stability	0			Physical hazards	0
	Physical and chemical properties	-			Personal protection	-

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

CAS No 483-84-1  
Molecular Weight 314.23  
Formula  $C_{10}H_6N_2O_8S$

Chemical Name	CAS No	Weight %	Oral LD50	Dermal LD50	Inhalation LC50
Flavianic acid hydrate	483-84-1	>98	-	-	-

## 4. FIRST AID MEASURES

### First Aid Measures

General advice	Immediate medical attention is required.
Eye contact	Immediate medical attention is required. Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Keep eye wide open while rinsing. Do not rub affected area.
Skin Contact	Immediate medical attention is required. Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes.
Inhalation	Remove to fresh air Call a physician or poison control center immediately If not breathing, give artificial respiration If breathing is difficult, give oxygen
Ingestion	Immediate medical attention is required. Do NOT induce vomiting. Drink plenty of water. Never give anything by mouth to an unconscious person. Remove from exposure, lie down. Clean mouth with water and drink afterwards plenty of water. Call a physician or poison control center immediately.
Self-protection of the first aider	Use personal protective equipment as required. Avoid contact with skin, eyes or clothing.

### Most important symptoms and effects, both acute and delayed

Symptoms No information available.

### Indication of any immediate medical attention and special treatment needed

Note to physicians Product is a corrosive material. Use of gastric lavage or emesis is contraindicated. Possible perforation of stomach or esophagus should be investigated. Do not give chemical antidotes. Asphyxia from glottal edema may occur. Marked decrease in blood pressure may occur with moist rales, frothy sputum, and high pulse pressure. Treat symptomatically.

## 5. FIRE-FIGHTING MEASURES

### Suitable Extinguishing Media

---

Suitable Extinguishing Media	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
Unsuitable Extinguishing Media	None.

### **Specific hazards arising from the chemical**

Specific hazards arising from the chemical	The product causes burns of eyes, skin and mucous membranes. Thermal decomposition can lead to release of irritating and toxic gases and vapors. In the event of fire and/or explosion do not breathe fumes.
Hazardous combustion products	Carbon oxides. Nitrogen oxides (NOx).

### **Explosion data**

Sensitivity to Mechanical Impact	No information available.
Sensitivity to Static Discharge	No information available.

### **Protective equipment and precautions for firefighters**

Protective equipment and precautions for firefighters	As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.
---	--

## **6. ACCIDENTAL RELEASE MEASURES**

---

### **Personal precautions, protective equipment and emergency procedures**

Personal precautions	Evacuate personnel to safe areas. Use personal protective equipment as required. Avoid contact with skin, eyes or clothing. Keep people away from and upwind of spill/leak.
----------------------	---

### **Environmental precautions**

Environmental precautions	Do not allow into any sewer, on the ground or into any body of water. Should not be released into the environment. Prevent further leakage or spillage if safe to do so. Prevent product from entering drains. See Section 12 for additional Ecological Information.
---------------------------	--

### **Methods and material for containment and cleaning up**

Methods for containment	Prevent further leakage or spillage if safe to do so.
Methods for cleaning up	Take up mechanically, placing in appropriate containers for disposal. Clean contaminated surface thoroughly. Use personal protective equipment as required. Cover powder spill with plastic sheet or tarp to minimize spreading and keep powder dry. Avoid creating dust. Dam up. After cleaning, flush away traces with water.

## **7. HANDLING AND STORAGE**

---

### **Precautions for safe handling**

Advice on safe handling	Use personal protective equipment as required. Avoid contact with skin, eyes or clothing. Ensure adequate ventilation, especially in confined areas. In case of insufficient ventilation, wear suitable respiratory equipment. Use only with adequate ventilation and in closed systems.
-------------------------	--

### **Conditions for safe storage, including any incompatibilities**

Storage Conditions	Keep container tightly closed in a dry and well-ventilated place. Keep out of the reach of children. Keep containers tightly closed in a dry, cool and well-ventilated place. Keep in properly labeled containers. Store at room temperature.
Incompatible materials	Incompatible with strong acids and bases. Incompatible with oxidizing agents.

## **8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

---

---

## Control parameters

### Exposure Guidelines

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

## Appropriate engineering controls

### Engineering Controls

Showers  
Eyewash stations  
Ventilation systems

## Individual protection measures, such as personal protective equipment

### Eye/face protection

Tight sealing safety goggles. Face protection shield.

### Skin and Body Protection

Wear protective gloves and protective clothing.

### Respiratory protection

If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

### General Hygiene Considerations

When using do not eat, drink or smoke. Wash contaminated clothing before reuse. Keep away from food, drink and animal feeding stuffs. Contaminated work clothing should not be allowed out of the workplace. Regular cleaning of equipment, work area and clothing is recommended. Avoid contact with skin, eyes or clothing. Take off all contaminated clothing and wash it before reuse. Wear suitable gloves and eye/face protection.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

### Physical State

Solid

### Appearance

crystalline or powder

### Odor

No information available

### Property

### Values

#### pH

No information available

#### Melting point/freezing point

148 °C

#### Boiling point

No information available

#### Flash point

No information available

#### Density

No information available

#### Evaporation rate

No information available

#### Upper flammability limits

No information available

#### Lower flammability limit

No information available

#### Vapor pressure

No information available

#### Vapor density

No information available

#### Specific gravity

No information available

#### Water solubility

No information available

#### Solubility in other solvents

No information available

#### Partition coefficient

No information available

#### Autoignition temperature

No information available

#### Decomposition temperature

No information available

#### Kinematic viscosity

No information available

#### Explosive properties

No information available

#### Oxidizing properties

No information available

---

## 10. STABILITY AND REACTIVITY

### Reactivity

Not applicable

### Chemical stability

Stable under recommended storage conditions.

### Possibility of Hazardous Reactions

None under normal processing.

### Hazardous polymerization

No information available.

### Conditions to avoid

Exposure to air or moisture over prolonged periods.

### Incompatible materials

Incompatible with strong acids and bases. Incompatible with oxidizing agents.

---

Hazardous Decomposition Products	Thermal decomposition can lead to release of irritating and toxic gases and vapors. Carbon oxides. Nitrogen oxides (NOx).
----------------------------------	---

## 11. TOXICOLOGICAL INFORMATION

---

### Information on likely routes of exposure

Inhalation	No data available.
Eye contact	No data available.
Skin Contact	No data available.
Ingestion	No data available.

### Information on toxicological effects

Symptoms	No information available.
----------	---------------------------

### Delayed and immediate effects as well as chronic effects from short and long-term exposure

Chronic Toxicity	Chronic exposure to corrosive fumes/gases may cause erosion of the teeth followed by jaw necrosis. Bronchial irritation with chronic cough and frequent attacks of pneumonia are common. Gastrointestinal disturbances may also be seen. Avoid repeated exposure. Possible risk of irreversible effects.
------------------	--

### Numerical measures of toxicity - Product Information

Unknown acute toxicity	No information available
------------------------	--------------------------

## 12. ECOLOGICAL INFORMATION

---

Ecotoxicity	May cause long lasting harmful effects to aquatic life
-------------	--

100% of the mixture consists of components(s) of unknown hazards to the aquatic environment.

Persistence and degradability	No information available.
Bioaccumulation	No information available.
Mobility	No information available.

## 13. DISPOSAL CONSIDERATIONS

---

Disposal of wastes	Disposal should be in accordance with applicable regional, national and local laws and regulations.
Contaminated packaging	Do not reuse container.
Other Information	Waste codes should be assigned by the user based on the application for which the product was used.
US EPA Waste Number	D002

## 14. TRANSPORT INFORMATION

---

### DOT

UN/ID no	UN1759
Hazard Class	8
Packing Group	II
Proper shipping name	Corrosive solid, n.o.s

---

Description UN1759, Corrosive solid, n.o.s (Flavianic acid hydrate), 8, II  
Emergency Response Guide Number 154

IMDG Not regulated

IATA Not regulated

## 15. REGULATORY INFORMATION

---

### International Inventories

**All of the components in the product are on the following Inventory lists**

TSCA (United States): Canada (DSL/NDSL) Europe (EINECS/ELINCS/NLP) South Korea (KECL):

Chemical Name	TSCA	DSL	NDSL	EINECS	ELINCS	ENCS	IECSC	KECL	PICCS	AICS
Flavianic acid hydrate	X	-	X	X	-	-	-	X	-	-

*X - Listed*

*TSCA - United States Toxic Substances Control Act Section 8(b) Inventory*

*DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List*

*EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances*

*ENCS - Japan Existing and New Chemical Substances*

*IECSC - China Inventory of Existing Chemical Substances*

*KECL - Korean Existing and Evaluated Chemical Substances*

*PICCS - Philippines Inventory of Chemicals and Chemical Substances*

### US Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

#### **SARA 311/312 Hazard Categories**

Acute health hazard	Yes
Chronic Health Hazard	Yes
Fire hazard	No
Sudden release of pressure hazard	No
Reactive hazard	No

#### **CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

### US State Regulations

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

#### **U.S. State Right-to-Know Regulations**

This product does not contain any substances regulated by state right-to-know regulations

## 16. OTHER INFORMATION

---

Revision note No information available

---

**Disclaimer**

The information provided in this Material Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.