

Preparation Date: 01/01/2019

Revision Date: N/A

Revision Number: N/A

## 1. IDENTIFICATION

### Product identifier

**Product code:** C6470  
**Product Name:** POTASSIUM IODIDE, GRANULAR, REAGENT, ACS

### Other means of identification

**Synonyms:** No information available  
**CAS #:** 7681-11-0  
**RTECS #** TT2975000  
**CI#:** Not available

### Recommended use of the chemical and restrictions on use

**Recommended use:** No information available.  
**Uses advised against** No information available

**Supplier:** **Dawn Scientific Inc**  
121 Liberty Street, Metuchen, NJ, 08840  
Tel : 732-902-6300 | Fax : 973-802-1005  
sales@dawnscientific.com | www.dawnscientific.com

**Emergency telephone number** Chemtrec 1-800-424-9300

## 2. HAZARDS IDENTIFICATION

### Classification

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Considered a dangerous substance or mixture according to the Globally Harmonized System (GHS)

Serious eye damage/eye irritation	Category 2B
Skin sensitization	Category 1B
Reproductive toxicity	Category 2
Specific target organ toxicity (single exposure)	Category 3
Specific target organ toxicity (repeated exposure)	Category 1

### Label elements

#### **Danger**

#### **Hazard statements**

Causes eye irritation  
May cause an allergic skin reaction  
Suspected of damaging fertility or the unborn child  
May cause respiratory irritation

Causes damage to organs through prolonged or repeated exposure



**Hazards not otherwise classified (HNOC)**

Not Applicable

**Other hazards**

Not available

**Precautionary Statements - Prevention**

Obtain special instructions before use  
Do not handle until all safety precautions have been read and understood  
Wash face, hands and any exposed skin thoroughly after handling  
Wear protective gloves/protective clothing/eye protection/face protection  
Contaminated work clothing must not be allowed out of the workplace  
Use only outdoors or in a well-ventilated area  
Do not breathe dust/fume/gas/mist/vapors/spray  
Do not eat, drink or smoke when using this product

**Precautionary Statements - Response**

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical attention.  
IF ON SKIN: Wash with plenty of water  
If skin irritation or rash occurs: Get medical attention  
Wash contaminated clothing before reuse  
IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER or physician if you feel unwell.

**Precautionary Statements - Storage**

Store locked up  
Store in a well-ventilated place. Keep container tightly closed

**Precautionary Statements - Disposal**

Dispose of contents and container to an approved waste disposal plant in accordance with local, regional, national and international regulations as applicable

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

Component	CAS No	Weight-%
Potassium Iodide	7681-11-0	100

### 4. FIRST AID MEASURES

**First aid measures**

**General Advice:**

National Capital Poison Center in the United States can provide assistance if you have a poison emergency and need to talk to a poison specialist. Call 1-800-222-1222.

**Skin Contact:**

Wash off immediately with soap and plenty of water removing all contaminated clothing and shoes. Get medical attention if irritation develops. Consult a physician if necessary.

<b>Eye Contact:</b>	Flush eyes with water for 15 minutes. Get medical attention if irritation occurs. If symptoms persist, call a physician.
<b>Inhalation:</b>	Move to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention.
<b>Ingestion:</b>	Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person. Consult a physician if necessary.

#### **Most important symptoms and effects, both acute and delayed**

<b>Symptoms</b>	<p>Causes eye irritation</p> <p>Eye contact may result in redness or pain</p> <p>May cause an allergic skin reaction</p> <p>May cause irritation of respiratory tract</p> <p>Causes digestive (gastrointestinal) tract irritation</p> <p>Ingestion may cause nausea, vomiting, and diarrhea</p> <p>May affect respiration</p> <p>Dyspnea (Shortness of breath and difficulty breathing)</p> <p>It may affect the thyroid</p> <p>Central nervous system effects</p> <p>May cause headache</p> <p>Drowsiness</p> <p>Somnolence</p> <p>May cause iodism</p>
-----------------	--

#### **Indication of any immediate medical attention and special treatment needed**

<b>Notes to Physician:</b>	Treat symptomatically.
----------------------------	------------------------

#### **Protection of first-aiders**

First-Aid Providers: Avoid exposure to blood or body fluids. Wear gloves and other necessary protective clothing. Dispose of contaminated clothing and equipment as bio-hazardous waste.

### **5. FIRE-FIGHTING MEASURES**

#### **Extinguishing Media**

<b>Suitable Extinguishing Media:</b>	The product is not flammable. If it is involved in a fire, extinguish the fire using an agent suitable for the type of surrounding fire.
--------------------------------------	--

<b>Unsuitable Extinguishing Media:</b>	No information available.
--	---------------------------

#### **Specific hazards arising from the chemical**

<b>Hazardous combustion products</b>	No information available.
--------------------------------------	---------------------------

<b>Specific hazards</b>	No information available.
-------------------------	---------------------------

#### **Special Protective Actions for Firefighters**

<b>Specific Methods:</b>	No information available
--------------------------	--------------------------

<b>Special Protective Equipment for Firefighters:</b>	As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear
---	---

## 6. ACCIDENTAL RELEASE MEASURES

### Personal precautions, protective equipment and emergency procedures

**Personal Precautions:** Ensure adequate ventilation. Keep people away from and upwind of spill/leak. Use personal protective equipment. Avoid contact with skin, eyes and clothing. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing. Avoid breathing dust.

**Environmental precautions** Prevent further leakage or spillage if safe to do so. Prevent product from entering drains. Prevent entry into waterways, sewers, basements or confined areas.

### Methods and material for containment and cleaning up

**Methods for containment** Stop leak if you can do it without risk. Cover with plastic sheet to prevent spreading.

**Methods for cleaning up** Sweep up and shovel into suitable containers for disposal. Clean contaminated surface thoroughly.

## 7. HANDLING AND STORAGE

### Precautions for safe handling

#### **Technical Measures/Precautions:**

Provide sufficient air exchange and/or exhaust in work rooms. Avoid dust formation. Keep away from incompatible materials.

#### **Safe Handling Advice:**

Wear personal protective equipment. Avoid contact with skin, eyes and clothing. Keep away from heat and sources of ignition. Avoid dust formation. Do not ingest. Do not breathe vapors/dust. Handle in accordance with good industrial hygiene and safety practice.

### Conditions for safe storage, including any incompatibilities

#### **Technical Measures/Storage Conditions:**

Keep container tightly closed in a dry and well-ventilated place. Store at room temperature in the original container. Store away from incompatible materials. Protect from moisture. Moisture sensitive. Protect from light. Sensitive to light. Store in light-resistant containers.

#### **Incompatible Materials:**

Alkali Metals  
Metals  
Organic materials  
Reducing agents  
Strong acids  
Oxidizing agents  
Mercurous chloride

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### Control parameters

### **National occupational exposure limits**

#### **United States**

Component	CAS No	OSHA	NIOSH	ACGIH	AIHA WEEL
Potassium Iodide	7681-11-0	None	None	0.01 ppm TWA	None

				inhalable fraction and vapor (for iodides)	
--	--	--	--	--	--

#### Canada

Component	CAS No	Canada - Alberta	Canada - British Columbia	Canada - Ontario	Canada - Quebec
Potassium Iodide	7681-11-0	None	None	None	None

#### Australia and Mexico

Component	CAS No	Australia	Mexico
Potassium Iodide	7681-11-0	None	None

#### Appropriate engineering controls

##### Engineering measures to reduce exposure:

Ensure adequate ventilation. Use process enclosures, local exhaust ventilation, or other engineering controls to keep airborne levels below recommended exposure limits. If user operations generate dust, fume or mist, use ventilation to keep exposure to airborne contaminants below the exposure limit.

#### Individual protection measures, such as personal protective equipment

##### Personal Protective Equipment

<b>Eye protection:</b>	Safety glasses with side-shields. or Goggles
<b>Skin and body protection:</b>	Long sleeved clothing Chemical resistant apron Gloves
<b>Respiratory protection:</b>	Effective dust mask. Use a dust respirator under conditions where exposure to the substance is apparent (e.g. generation of high concentration of dust (dust clouds) , inadequate ventilation, development of respiratory tract irritation), and engineering controls are not feasible. Be sure to use an approved/certified respirator or equivalent.
<b>Hygiene measures:</b>	Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product When using, do not eat, drink or smoke.

### 9. PHYSICAL AND CHEMICAL PROPERTIES

<b>Physical state:</b> Solid	<b>Appearance:</b> Crystals. Crystalline. Granular. Granular powder.	<b>Color:</b> White.
<b>Odor:</b> Odorless.	<b>Taste</b> Bitter. Saline.	<b>Formula</b> KI
<b>Molecular/Formula weight (g/mole):</b> 166.00	<b>Flammability (solid, gas)</b> no data available	<b>Flashpoint (°C/°F):</b> No information available
<b>Flash Point Tested according to:</b> Not available	<b>Autoignition Temperature (°C/°F):</b> No information available	<b>Lower Explosion Limit (%):</b> No information available
<b>Upper Explosion Limit (%):</b>		

No information available	<b>Melting point/range(°C/°F):</b> 681°C/1257.8°F	<b>Decomposition temperature(°C/°F):</b> No information available
<b>Boiling point/range(°C/°F):</b> 1330°C/2426°F	<b>Bulk density:</b> No information available	<b>Density (g/cm3):</b> No information available
<b>Specific gravity:</b> 3.1	<b>pH</b> No information available	<b>Vapor pressure @ 20°C (kPa):</b> No information available
<b>Evaporation rate:</b> No information available	<b>Vapor density:</b> No information available	<b>VOC content (g/L):</b> No information available
<b>Odor threshold (ppm):</b> No information available	<b>Partition coefficient (n-octanol/water):</b> No information available	<b>Viscosity:</b> No information available
<b>Miscibility:</b> No information available	<b>Solubility:</b> Easily soluble in cold water Easily soluble in water Soluble in Methanol Partially soluble in acetone	

## 10. STABILITY AND REACTIVITY

### Reactivity

Reactive with acids

Reactive with metals

Reactive with oxidizing agents

Reacts with reducing agents

Reactive with organic materials

Mercurous chloride in the presence of an excess of potassium iodide produces metallic mercury and mercuric iodide, the latter forming the soluble double salt, potassium mercuric iodide

Reacts violently with strong oxidizers (bromotrifluorides, chlorotrifluorides, fluorine perchlorate), metallic salts. Attacks metals in moist environments. Also incompatible with salts of alkaloids, chloral hydrate, calomel (mercurous chloride), potassium chlorate, tartaric and other acids, oxidants, diazonium salts, charcoal, ozone, strong reducers, alkali metals, metals (brass, aluminum magnesium, zinc, cadmium, copper, tin, nickel, steel)

### Chemical stability

**Stability:** Stable under recommended storage conditions.

**Possibility of Hazardous Reactions:** Hazardous polymerization does not occur

**Conditions to avoid:** Heat. Ignition sources. Incompatible materials. Exposure to light. Exposure to moisture. Becomes yellow on exposure to air. Light and moisture accelerate the decomposition.

**Incompatible Materials:** Alkali Metals  
Metals  
Organic materials  
Reducing agents  
Strong acids  
Oxidizing agents  
Mercurous chloride

**Hazardous decomposition products:** Hydrogen iodide. Oxides of potassium. Iodine.

### Other Information

**Corrosivity:** No information available

**Special Remarks on Corrosivity:** No information available

## 11. TOXICOLOGICAL INFORMATION

### Information on likely routes of exposure

**Principal Routes of Exposure:**

Ingestion. Inhalation.

### Acute Toxicity

### Component Information

Potassium Iodide	
CAS No	7681-11-0

**LD50/oral/rat** = No information available

**LD50/oral/mouse** = 1862 mg/kg, LDLo

**LD50/dermal/rabbit** = No information available

**LD50/dermal/rat** = No information available

**LC50/inhalation/rat** = No information available

**LC50/inhalation/mouse** = No information available

**Other LD50 or LC50 information** = 916 mg/kg, oral, rabbit, LDLo

### Product Information

**LD50/oral/rat** =

**Value - Acute Toxicity** = No information available

**LD50/oral/mouse** =

**Value - Acute Tox** = No information available

**LD50/dermal/rabbit**

**Value - Acute Toxicity** = No information available

**LD50/dermal/rat**

**VALUE - Acute Tox** = No information available

**LC50/inhalation/rat**

**VALUE-Vapor** = No information available

**VALUE-Gas** = No information available

**VALUE-Dust/Mist** = No information available

**LC50/Inhalation/mouse**

**VALUE-Vapor** = No information available

**VALUE - Gas** = No information available

**VALUE - Dust/Mist** = No information available

### Symptoms

**Skin Contact:** May cause skin irritation. It may cause mild skin irritation. May cause allergic contact dermatitis.

**Eye Contact:** Causes eye irritation. Mild eye irritation.

**Inhalation** May cause respiratory tract and mucous membrane irritation and a productive cough. May cause pulmonary edema and inflammation of the tonsils.

**Ingestion** Causes gastrointestinal tract irritation with nausea, vomiting and diarrhea. May affect behavior(somnolence, muscle weakness), respiration (dyspnea). Serum-sickness type of hypersensitivity such as fever, arthralgia, lymph node enlargement, and eosinophilia may appear. Thrombotic thrombocytopenic purpura, and fatal periarteritis nodosa attributed to hypersensitivity to iodide has been described.

**Aspiration hazard** No information available.

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

**Chronic Toxicity** Can lead to iodism characterized by salivation, metallic taste, sore mouth, nausea, vomiting, epigastric pain, diarrhea, nasal discharge, sneezing, conjunctivitis, fever, headache, laryngitis, bronchitis, stomatitis, parotitis, anemia, and skin rashes. Chronic ingestion may also affect metabolism (anorexia), and thyroid gland (hypothyroidism, goiter). Furthermore, chronic ingestion of iodides (in animals) during pregnancy has resulted in fetal deaths, severe goiter and cretinoid appearance of the newborn. Prolonged or repeated skin contact may cause allergic reaction.

**Sensitization:** May cause sensitization by skin contact.

**Mutagenic Effects:** May affect genetic material based on animal test data  
No human data found

**Carcinogenic effects:** Not considered carcinogenic.

Component	CAS No	IARC	ACGIH - Carcinogens	NTP	OSHA HCS - Carcinogens	Australia - Notifiable Carcinogenic Substances	Australia - Prohibited Carcinogenic Substances
Potassium Iodide	7681-11-0	Not listed	A4 - Not Classifiable as a Human Carcinogen (listed as Iodides)	Not listed	Not listed	Not listed	Not listed

ACGIH (American Conference of Governmental Industrial Hygienists)

IARC (International Agency for Research on Cancer)

NTP (National Toxicology Program)

OSHA (Occupational Safety and Health Administration of the US Department of Labor)

**Reproductive toxicity** Suspected of damaging fertility or the unborn child

**Reproductive Effects:** No information available

**Developmental Effects:** No information available

**Teratogenic Effects:** No information available

**Specific Target Organ Toxicity**

**STOT - single exposure** respiratory system.

**STOT - repeated exposure** Causes damage to organs through prolonged or repeated exposure.

**Target Organs:** Thyroid.



## 12. ECOLOGICAL INFORMATION

### Ecotoxicity

<b>Ecotoxicity effects:</b>	Aquatic environment.
<i>Potassium Iodide - 7681-11-0</i>	
<b>Fish</b>	LC50: >100mg/L (96h, Danio rerio)
<b>Persistence and degradability:</b>	No information available
<b>Bioaccumulative potential:</b>	No information available.
<b>Mobility in soil</b>	No information available
<b>Other adverse effects</b>	No information available.

## 13. DISPOSAL CONSIDERATIONS

### Disposal Methods

#### Waste from residues / unused products:

Waste must be disposed of in accordance with Federal, State and Local regulation.

#### Contaminated packaging:

Empty containers should be taken for local recycling, recovery or waste disposal

Component	CAS No	RCRA - F Series Wastes	RCRA - K Series Wastes	RCRA - P Series Wastes	RCRA - U Series Wastes
Potassium Iodide	7681-11-0	None	None	None	None

## 14. TRANSPORT INFORMATION

### DOT

<b>UN-No:</b>	Not Regulated
<b>Proper Shipping Name:</b>	No information available
<b>Hazard Class</b>	No information available
<b>Subsidiary Class</b>	No information available
<b>Packing group:</b>	No information available
<b>Emergency Response Guide Number</b>	No information available
<b>Marine Pollutant</b>	No data available
<b>DOT RQ (lbs):</b>	No information available
<b>Special Provisions</b>	No Information available
<b>Symbol(s):</b>	No information available
<b>Description:</b>	No information available

### TDG (Canada)

<b>UN-No:</b>	Not Regulated
<b>Proper Shipping Name:</b>	No information available
<b>Hazard Class</b>	No information available
<b>Subsidiary Risk:</b>	No information available
<b>Packing Group:</b>	No information available
<b>Marine Pollutant</b>	No Information available
<b>Description:</b>	No information available

**ADR**

UN Number Not regulated  
 Proper Shipping Name: No information available  
 Transport hazard class(es) No information available  
 Packing group No information available  
 Subsidiary Risk: No information available

**IMDG**

UN-No: Not Regulated  
 Proper Shipping Name: No information available  
 Hazard Class: No information available  
 Subsidiary Risk: No information available  
 Packing Group: No information available  
 Marine Pollutant No information available

**RID**

UN Number Not Regulated  
 Proper Shipping Name: No information available  
 Transport hazard class(es) No information available  
 Subsidiary Risk: No information available  
 Packing group No information available

**ICAO (air)**

UN-No: Not Regulated  
 Proper Shipping Name: No information available  
 Hazard Class No information available  
 Subsidiary Risk: No information available  
 Packing Group: No information available

**IATA**

UN Number Not Regulated  
 Proper Shipping Name: No information available  
 Transport hazard class(es) No information available  
 Subsidiary Risk: No information available  
 Packing group No information available  
 Precautionary Statements - No information available  
 Response  
 Special Provisions No information available

<b>15. REGULATORY INFORMATION</b>
-----------------------------------

**International Inventories**

Component	CAS No	U.S. TSCA	KOREA KECL	Philippines (PICCS)	Japan ENCS	China IECSC	Australia (AICS)	EINECS-No.
Potassium Iodide	7681-11-0	PresentACTIVE	Present KE-29149	Present	Present (1)-439	Present	Present	Present 231-659-4

**U.S. Regulations***Potassium Iodide*

FDA - Food Additives Generally Recognized as Safe (GRAS): 21 CFR 184.1634

FDA - Direct Food Additives 21 CFR 172.375

FDA - 21 CFR - Total Food Additives 172.375, 178.1010, 184.1634

- List Sourced from EAFUS

California Prop. 65: Safe Drinking Water and Toxic Enforcement Act of 1986.

**Chemicals Known to the State of California to Cause Cancer:**

This product does not contain a chemical requiring a warning under California Prop. 65. (See table below)

**Chemicals Known to the State of California to Cause Reproductive Toxicity:**

This product does not contain a chemical requiring a warning under California Prop. 65. (See table below)

Component	CAS No	Carcinogen	Developmental Toxicity	Male Reproductive Toxicity	Female Reproductive Toxicity:
Potassium Iodide	7681-11-0	Not Listed	Not Listed	Not Listed	Not Listed

**CERCLA/SARA**

Component	CAS No	CERCLA - Hazardous Substances and their Reportable Quantities	Section 302 Extremely Hazardous Substances and TPQs	Section 302 Extremely Hazardous Substances and RQs	Section 313 - Chemical Category	Section 313 - Reporting de minimis
Potassium Iodide	7681-11-0	None	None	None	None	None

**U.S. TSCA**

Component	CAS No	TSCA Section 5(a)2 - Chemicals With Significant New Use Rules (SNURS)	TSCA 8(d) -Health and Safety Reporting
Potassium Iodide	7681-11-0	Not Applicable	Not Applicable

**Canada****WHMIS 2015 - GHS Classifications**

WHMIS 2015 Hazard Classification Information:

Component  
Potassium Iodide  
7681-11-0 ( 100 )

WHMIS 2015 Hazard Classification  
Specific target organ toxicity - Repeated exposure - Category 1:  
H372 Causes damage to organs through prolonged or repeated exposure.

**Canada Hazardous Products Regulation** This product has been classified according to the hazard criteria of the HPR (Hazardous Products Regulation) and the SDS contains all of the information required by the HPR

**DSL/NDSL**

Component	CAS No	Canada (DSL)	Canada (NDSL)
Potassium Iodide	7681-11-0	Present	Not Listed

Component	CAS No	CEPA Schedule I - Toxic Substances
Potassium Iodide	7681-11-0	Not listed
Component	CAS No	CEPA - 2010 Greenhouse Gases Subject to Mandatory Reporting
Potassium Iodide	7681-11-0	Not listed

**EU Classification****EU GHS - SV - CLP 1272/2008**

Component	CAS No	EU GHS - SV - CLP (1272/2008)
-----------	--------	-------------------------------

Potassium Iodide	7681-11-0	No information
------------------	-----------	----------------

EU - CLP (1272/2008)

**R-phrase(s)**

Not determined  
none

**S -phrase(s)**

none

Component	CAS No	Classification	Concentration Limits:	Safety Phrases
Potassium Iodide	7681-11-0		No information	

The product is classified in accordance with Annex VI to Directive 67/548/EEC

**Indication of danger:**

not determined

## 16. OTHER INFORMATION

**Preparation Date:** 01/01/2019  
**Revision date** N/A  
**Prepared by:** -

**Disclaimer:**

All chemicals may pose unknown hazards and should be used with caution. This Safety Data Sheet (SDS) applies only to the material as packaged. If this product is combined with other materials, deteriorates, or becomes contaminated, it may pose hazards not mentioned in this SDS. The physical properties reported in this SDS are obtained from the literature and do not constitute product specifications. Information contained herein does not constitute a warranty, whether expressed or implied, as to the safety, merchantability or fitness of the goods for a particular purpose. Dawn Scientific Inc Chemicals & Laboratory Products, assumes no responsibility for results obtained or for incidental or consequential damages, including lost profits, arising from the use of these data. No warranty against infringement of any patent, copyright or trademark is made or implied. It shall be the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use. While this SDS is based on technical data judged to be reliable, Dawn Scientific Inc assumes no responsibility for the completeness or accuracy of the information contained herein.

**End of Safety Data Sheet**