



Wide-Tip Volumetric Pipets

Reusable Glass, To Deliver

SIBATA 2040-Series Wide-Tip Volumetric Pipets meet ASTM E-969 Class B Specifications for Tolerance.

The wide-tip opening facilitates rapid fill and drain of viscous liquids, slurries, and suspensions; an excellent product for water treatment facilities, winemakers, and essential oil testing.

These pipets have sharp and permanent graduation lines and are calibrated to deliver. Feature permanent amber markings.

Sizes 1, 2, & 5mL do not have drawn tips, while sizes 10mL through 100mL feature standard drawn tips.

Specifications

| Cat. No. | Capacity (mL) | Tolerance (\pm mL) | Approx. Length (mm) | Approx. Tip I.D. (mm) |
|--------------|---------------|-----------------------|---------------------|-----------------------|
| 2040-1-1.9MM | 1.0 | 0.01 | 320 | 1.9 |
| 2040-2-2.4MM | 2.0 | 0.01 | 330 | 2.4 |
| 2040-5-3MM | 5.0 | 0.02 | 346 | 3.0 |
| 2040-10-3MM | 10.0 | 0.02 | 440 | 3.0 |
| 2040-20-3MM | 20.0 | 0.03 | 478 | 3.0 |
| 2040-25-3MM | 25.0 | 0.03 | 475 | 3.0 |
| 2040-50-3MM | 50.0 | 0.05 | 538 | 3.0 |
| 2040-100-3MM | 100.0 | 0.10 | 605 | 3.0 |

Packaging

| Cat. No. | Pack Qty. | Case Qty. |
|--------------|-----------|-----------|
| 2040-1-1.9MM | 12 | 72 |
| 2040-2-2.4MM | 12 | 72 |
| 2040-5-3MM | 12 | 72 |
| 2040-10-3MM | 12 | 72 |
| 2040-20-3MM | 6 | 36 |
| 2040-25-3MM | 6 | 36 |
| 2040-50-3MM | 6 | 36 |
| 2040-100-3MM | 6 | 36 |

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



Available
at
Comparable
Price



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCI

FOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM