



DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



Graduated Cylinders

Student, Plastic Base, Double Metric Scale, Graduated, To Contain

SIBATA 2355-Series Graduated Cylinders meet ASTM E-1272 Class B, Style I Specifications for Classification and Tolerances.

These cylinders are made of two parts: a strong, accurate glass graduate and a detachable, sturdy plastic hexagonal base which absorbs impacts and reduces breakage. Feature hexagonal bases for added stability. All cylinders come with protective, polyethylene, hexagonal bumper guards.

All markings are in durable, white enamel.

Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Tolerance (\pm mL)	Approx. O.D. x Height (mm)
2355-10	10	0.4 to 10.0	0.1	0.10	13 x 177
2355-25	25	1.4 to 25.0	0.2	0.30	16 x 210
2355-50	50	2.0 to 50.0	1.0	0.40	24 x 212
2355-100	100	5.0 to 100.0	1.0	0.60	30 x 232

Packaging

Cat. No.	Pack Qty.	Case Qty.
2355-10	6	24
2355-25	6	24
2355-50	6	24
2355-100	6	24



Available at Comparable Price

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCI

FOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>