

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER



## Volumetric Flasks

Class A, Low-Actinic Amber,  $\Phi$  Stopper, To Contain

SIBATA 2307A-Series Low-Actinic Amber Volumetric Flasks meet ASTM E-288 Class A Specifications & USP Standards for Volumetric Glassware.

These flasks feature machine-blown bodies and heavy-beaded, heavy-tubing necks, tooled for  $\Phi$  stoppers. The graduation line is sharp and permanent; large, white, block-letters allow the labeling to be easily read. Each flask also has a large marking square. The low-actinic amber glass provides protection for materials sensitive to light.

All sizes are supplied with a glass  $\Phi$  stopper.


### Specifications

Cat. No.	Capacity (mL)	Tolerance ( $\pm$ mL)	Stopper No.	Approx. O.D. x Height (mm)
2307A-10	10	0.02	9	28 x 95
2307A-20	20	0.03	9	36 x 115
2307A-25	25	0.03	9	39 x 122
2307A-50	50	0.05	13	49 x 150
2307A-100	100	0.08	13	60 x 180
2307A-200	200	0.10	16	75 x 220
2307A-250	250	0.12	16	78 x 239
2307A-500	500	0.15	19	100 x 280
2307A-1000	1000	0.30	22	127 x 325

### Packaging

Cat No.	Pack Qty.	Case Qty.
2307A-10	6	72
2307A-20	6	72
2307A-25	6	72
2307A-50	4	48
2307A-100	4	48
2307A-200	4	24
2307A-250	4	24
2307A-500	2	12
2307A-1000	2	8

Available  
at  
Comparable  
Price

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

**1-800-DAWN-SCI**

FOR INQUIRY

**+1 732-902-6300**

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."