

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER



Graduated Cylinders

Class B, Double Metric Scale, Graduated, To Contain

SIBATA 2351-Series Graduated Cylinders meet ASTM E-1272 Class B, Style I Specifications.

Feature hexagonal bases for added stability. All cylinders, with the exception of the 2000mL size, come with protective, polyethylene, hexagonal bumper guards.

All markings are in durable, white enamel.

Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Tolerance (\pm mL)	Approx. O.D. x Height (mm)
2351-10	10	0.7 to 10.0	0.1	0.10	8 x 182
2351-25	25	1.4 to 25.0	0.2	0.30	17 x 210
2351-50	50	3.5 to 50.0	0.5	0.40	24 x 218
2351-100	100	5.0 to 100.0	1.0	0.60	30 x 236
2351-250	250	10.0 to 250.0	2.0	1.40	41 x 301
2351-500	500	25.0 to 500.0	5.0	2.60	53 x 362
2351-1000	1000	50.0 to 1000.0	10.0	5.00	70 x 422
2351-2000	2000	100.0 to 2000.0	20.0	10.00	90 x 478

Packaging

Cat. No.	Pack Qty.	Case Qty.
2351-10	4	48
2351-25	4	48
2351-50	4	48
2351-100	4	48
2351-250	4	24
2351-500	2	12
2351-1000	2	8

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

VEE GEE®

Available
at
Competable
Price

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCIFOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM<https://www.dawnscientific.com>