



DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



Filtration Filter Holders

Glass Microanalysis, 47mm

SIBATA 6168-Series Filter Holders are offered with a wide variety of components for different applications. Available with fritted glass support base for general and biological applications, or with removable stainless steel screen support base and PTFE gasket for ultraclean filtrate or proteinaceous solutions. The borosilicate glass funnels are available in the following capacity sizes:

Specifications

- Filter Diameter: 47mm
- Prefilter Diameter: 35mm
- Effective Filtration Area: 9.6cm²
- Maximum Capacity: 100mL Funnel: 150mL, 300mL Funnel: 350mL, 500mL Funnel: 550mL, 1000mL Funnel: 1100mL
- Graduated Volume: 100mL Funnel: 100mL, 300mL Funnel: 300mL, 500mL Funnel: 500mL, 1000mL Funnel: 1000mL
- Connection: Silicone Stopper (No. 8) or for 2000mL Filter Flask (No.9)
- Sterilization: Autoclavable (at 121°C / 250°F for 20 minutes)

- 100mL
- 300mL
- 500mL
- 1000mL



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCI

FOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."