



DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



Erlenmeyer Flasks

Wide Mouth, Low Form, Heavy-Duty Rim, Graduated

VEE GEE 20232-Series Erlenmeyer Flasks meet ASTM E-1404 Type I, Class 1 Specifications for Classification, Design, Capacity & Dimensions, and Markings.

These Erlenmeyer flasks feature heavy-duty rims to reduce chipping, white, permanent graduations, and large marking squares. Made from high-quality borosilicate tubing with a low coefficient of expansion. Borosilicate tubing also provides uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength.

Accuracy: $\pm 5\%$ of full capacity.

Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Rubber Stopper No.	Approx. O.D. Neck (mm)
20232-125	125	50 - 125	25	68 x 117	6	34
20232-250	250	100 - 250	25	84 x 138	8	50
20232-500	500	200 - 500	50	105 x 178	10	50
20232-1000	1000	400 - 1000	50	128 x 220	11	50
20232-2000	2000	600 - 1800	100	164 x 280	13	72

Packaging

Cat. No.	Pack Qty.	Case Qty.
20232-125	12	48
20232-250	12	48
20232-500	6	36
20232-1000	6	24
20232-2000	1	8



Available at Comparable Price

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

1-800-DAWN-SCI

FOR INQUIRY

+1 732-902-6300

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>