

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

VEE GEE®

SIBATA Class A Volumetric Flask Set

Sibata Class A Volumetric Flasks in a set of Four**One Each:**

- 10mL
- 25mL
- 50mL
- 100mL

SIBATA 2306A-Series Volumetric Flasks meet ASTM E-237 Class A, Style I Specifications. Feature machine-blown bodies and heavy-beaded, heavy-tubing necks, tooled for glass stoppers. The graduation line is sharp and precise. Large, white, block-letters allow the labeling to read. All flasks supplied with a glass stopper.



All markings in durable white enamel.

Specifications

Cat. No.	Capacity (mL)	Tolerance (± mL)	Stopper No.	Approx. O.D. x Height (mm)	Pack Qty.
202GK-004	10	0.02	9	26 x 96	1/PK/4
	25	0.03	9	37 x 123	
	50	0.05	13	49 x 152	
	100	0.08	13	58 x 181	

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

New!**Available
at
Comparable
Price****DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCIFOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM<https://www.dawnscientific.com>