

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

VEE GEE®

SIBATA Class A Graduated Cylinder Set

Sibata Class A Graduated Cylinders in a set of Four

One Each:

- 10mL
- 25mL
- 50mL
- 100mL



SIBATA 2351A-Series Graduated Cylinders meet ASTM E-1272 Class A, Style I Specifications & USP Standards for Volumetric Glassware. Feature hexagonal bases for added stability. All cylinders, with the exception of the 2000 ml size, feature protective, polyethylene, hexagonal bumper guards.

All markings are durable, white enamel.



Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Pack Qty.
202GK-006	10	0.7 - 10.0	0.1	8 x 182	1/PK/4
	25	1.4 - 25.0	0.2	17 x 210	
	50	3.5 - 50.0	0.5	24 x 218	
	100	5.0 - 100.0	1.0	30 x 236	

Available
at
Compitable
Price

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCIFOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM<https://www.dawnscientific.com>