

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

VEE GEE®

SIBATA Student Grade Graduated Cylinder Set

Sibata Student Grade Graduated Cylinders in a set of Four

One Each:

- 10mL
- 25mL
- 50mL
- 100mL

SIBATA 2355-Series Student Graduated Cylinders meet E-1272 Class B, Style I Specifications for Classification Tolerances. These cylinders are made of two parts: accurate glass graduate and a detachable, sturdy hexagonal base which absorbs impacts and reduces breakage. All cylinders feature protective hexagonal bumper guards



All markings are durable, white enamel.

Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Pack Qty.
202GK-005	10	0.4 - 10.0	0.1	13 x 177	1/PK/4
	25	1.4 - 25.0	0.2	16 x 210	
	50	2.0 - 50.0	1.0	24 x 212	
	100	5.0 - 100.0	1.0	30 x 232	

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

New!

Available
at
Competable
Price

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCIFOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM<https://www.dawnscientific.com>