

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

VEE GEE®

Filtering Flasks

Heavy Wall, Tubulation, Graduated

VEE GEE 20074-Series Filtering Flasks meet ASTM E-1406 Type II, Class 1 Specifications for Classification, Design, Capacity & Dimensions, and Markings.

These filtering flasks feature heavy-duty rims to reduce chipping, white, permanent graduations, and large marking squares. Made from high-quality borosilicate tubing with a low coefficient of expansion. Borosilicate tubing also provides uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength.

Tubulation outside diameter is 10mm. Accuracy: $\pm 5\%$ of full capacity.

Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Rubber Stopper No.
20074-250	250	75 - 250	25	86 X 160	6
20074-500	500	150 - 500	50	107 X 190	7
20074-1000	1000	400 - 1000	50	135 X 242	8
20074-2000	2000	600 - 2000	100	165 X 288	9

Packaging

Cat. No.	Pack Qty.	Case Qty.
20074-250	6	18
20074-500	6	18
20074-1000	6	12
20074-2000	1	8

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

1-800-DAWN-SCI

FOR INQUIRY

+1 732-902-6300

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>