



VEE GEE®

Burets

Class A, w/ PTFE Stopcocks

VEE GEE 20114-Series Burets meet ASTM E-287 Class A Specifications.

Made from precision, uniform Schott DURAN® true-bore tubing to ensure the stipulated capacity tolerances. These burets feature durable, permanent markings, fine, sharp lines, and large, easy-to-read numbers. Supplied with PTFE stopcocks to eliminate freezing and lubricant contamination and to provide a precision fit within the barrel.

All sizes are supplied with a 2mm stopcock bore plug.

Specifications

Cat. No.	Capacity (mL)	Graduation Interval (mL)	Tolerance (\pm mL)	Stopcock Bore (mm)	Approx. O.D. x Height (mm)
20114-10	10	0.05	0.02	2	9 x 510
20114-25	25	0.10	0.03	2	12 x 550
20114-50	50	0.10	0.05	2	13 x 758
20114-100	100	0.20	0.10	2	17 x 775

Packaging

Cat. No.	Pack Qty.	Case Qty.
20114-10	1	24
20114-25	1	24
20114-50	1	24
20114-100	1	24



"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

1-800-DAWN-SCI

FOR INQUIRY

+1 732-902-6300

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>