

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

VEE GEE®

Beakers

Griffin, Low-Form, Double Scale, Graduated

VEE GEE 20229-Series Beakers meet ASTM E-960 Type I Specifications for Classification, Design, Capacity & Dimensions, and Markings.

These beakers feature easy-pour spouts, white, permanent graduations, and large marking squares. Made from high-quality borosilicate glass with a low coefficient of expansion. Manufactured with uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength.

Cat. No.'s 20229-20 through 20229-3000 have double-graduated metric scales. Accuracy: $\pm 5\%$ of full capacity.

Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)
20229-10	10	--	--	25 x 35
20229-20	20	5 - 15	5	32 x 45
20229-30	30	5 - 25	5	34 x 55
20229-50	50	10 - 45	5	41 x 60
20229-100	100	20 - 80	10	50 x 75
20229-150	150	20 - 140	10	57 x 90
20229-250	250	25 - 200	25	68 x 100
20229-400	400	25 - 325	25	78 x 112
20229-600	600	50 - 500	50	90 x 125
20229-1000	1000	50 - 1000	50	108 x 150
20229-2000	2000	200 - 1800	100	130 x 188
20229-3000	3000	250 - 2500	125	149 x 212

Packaging

Cat. No.	Pack Qty.	Case Qty.
20229-10	12	48
20229-20	12	48
20229-30	12	48
20229-50	12	48
20229-100	12	48
20229-150	12	48
20229-250	12	48
20229-400	12	48
20229-600	6	36
20229-1000	6	24
20229-2000	4	8
20229-3000	1	6

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."



Available
at
Competable
Price

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

1-800-DAWN-SCI

FOR INQUIRY

+1 732-902-6300

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>