



VEE GEE®

Separatory Funnels

Squibb, Pear-Shaped, F Stopper, PTFE Stopcock

VEE GEE 20149-Series Separatory Funnels meet ASTM E-1096 Type 4 Specifications for Classification, Design, and Markings.

These separatory funnels feature smooth tooled tips, large, white block letters, and large marking squares. Supplied with PTFE stopcocks to eliminate freezing and lubricant contamination and to provide a precision fit within the barrel. The stopcock bores, body openings, and stems are carefully aligned for even liquid flow with fast cut-off. Made from Schott DURAN® borosilicate tubing with a low coefficient of expansion. Borosilicate tubing also provides uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength.

All sizes are supplied with a polyethylene F stopper.

Specifications

| Cat. No. | Capacity (mL) | F Stopper No. | Stopcock Plug Size | Approx. O.D. x Stem Length (mm) | Approx. Total Length (mm) |
|------------|---------------|------------------------|--------------------|---------------------------------|---------------------------|
| 20149-250 | 250 | 22 | 2 | 10 x 65 | 345 |
| 20149-500 | 500 | 27 | 4 | 10 x 65 | 392 |
| 20149-1000 | 1000 | 27 | 4 | 12 x 70 | 400 |
| 20149-2000 | 2000 | 38 | 6 | 12 x 75 | 570 |

Packaging

| Cat. No. | Pack Qty. | Case Qty. |
|------------|-----------|-----------|
| 20149-250 | 1 | 24 |
| 20149-500 | 1 | 24 |
| 20149-1000 | 1 | 12 |
| 20149-2000 | 1 | 9 |



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

1-800-DAWN-SCI

FOR INQUIRY

+1 732-902-6300

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>

DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



Available
at
Competable
Price