



**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER



# Erlenmeyer Set

Narrow Mouth, Low Form, Heavy-Duty Rim, Graduated Erlenmeyer Flasks available in a Three Piece Set.

One Each:

- 125mL
- 250mL
- 500mL

VEE GEE 20231-Series Erlenmeyer Flasks meet ASTM E-1404 Type I, Class 1 Specifications for Classification, Design, Capacity & Dimensions, and Markings.

These Erlenmeyer flasks feature heavy-duty rims to reduce chipping, white, permanent graduations, and large marking squares. Made from high-quality borosilicate tubing with a low coefficient of expansion. Borosilicate tubing also provides uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength.

Accuracy: ±5% of full capacity.



### Specifications

| Cat. No.  | Capacity (mL) | Graduation Range (mL) | Graduation Interval (mL) | Approx. O.D. x Height (mm) | Rubber Stopper No. | Approx. O.D. Neck (mm) | Pack Qty. |
|-----------|---------------|-----------------------|--------------------------|----------------------------|--------------------|------------------------|-----------|
| 202GK-003 | 125           | 50 - 125              | 25                       | 68 x 117                   | 5                  | 22                     | 1/PK/3    |
|           | 250           | 100 - 250             | 25                       | 84 x 138                   | 6                  | 34                     |           |
|           | 500           | 200 - 500             | 50                       | 105 x 178                  | 7                  | 34                     |           |

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

Available at Comparable Price



**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

**1-800-DAWN-SCI**

FOR INQUIRY

**+1 732-902-6300**

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>