



DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER

VEE GEE®

Beaker Set

VEE GEE, Low-Form Griffin Beakers available in a Five Piece Set.

One Each:

- 50mL
- 100mL
- 250mL
- 600mL
- 1000mL

VEE GEE 20229-Series Beakers meet ASTM E-960 Type I Specifications for Classification, Design, Capacity & Dimensions, and Markings.

These beakers feature easy-pour spouts, white, permanent graduations, and large marking squares. Made from high quality borosilicate glass with a low coefficient of expansion. Manufactured with uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength.

Accuracy: ±5% of full capacity.

Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Pack Qty.
202GK-002	50	10 - 45	5	41 x 60	1/PK/5
	100	20 - 80	10	50 x 75	
	250	25 - 200	25	68 x 100	
	600	50 - 500	50	90 x 125	
	1000	50 - 950	50	108 x 150	

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

1-800-DAWN-SCI

FOR INQUIRY

+1 732-902-6300

SALES@DAWNSCIENTIFIC.COM

<https://www.dawnscientific.com>