

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

**VEE GEE®**

# Media Bottles

VEE GEE Media Bottles are commonly used to contain chemicals in liquid or powder form in laboratories. The bottles are made from borosilicate glass in compliance with ASTM 438 Type 1 and are available in clear or amber (low actinic) for light sensitive chemical compounds.

Five sizes available in both clear and amber, VEE GEE media bottles have threaded on Blue Plastic caps and pouring rings that are autoclavable.

- Clear and Amber
- Made from Borosilicate Glass
- Threaded Blue Plastic Caps
- Autoclavable Caps & Pouring Ring
- Clear White Marking

## Media Bottles, Clear

Cat. No.	Capacity (mL)	Pack Qty.
20274-100	100	10
20274-250	250	10
20274-500	500	10
20274-1000	1000	10
20274-2000	2000	4

## Media Bottles, Amber

Cat. No.	Capacity (mL)	Pack Qty.
20275-100	100	10
20275-250	250	10
20275-500	500	10
20275-1000	1000	10
20275-2000	2000	4

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."



Available  
at  
Comparable  
Price

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER  
**1-800-DAWN-SCI**FOR INQUIRY  
**+1 732-902-6300**  
SALES@DAWNSCIENTIFIC.COM<https://www.dawnscientific.com>