

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

SIBATA SCIENTIFIC TECHNOLOGY LTD.



## Serological, Wide-Tip Pipets

### Reusable Glass, For Measuring

SIBATA wide-tip-series serological pipets for biochemical and other applications. Markings and color coding are permanently enameled to the glass. SIBATA pipets are meticulously manufactured for tolerances and accuracy. Compared to plastic, reusable glass pipets are autoclavable and sustainable with less waste. Glass is a superior material for measuring samples with a wide variety of chemical compounds.

### Three styles available:

- Wide-Tip, Two Lengths: 205mm-350mm
- Large Hole Opening, Short, 205mm Length
- Graduated to the Tip, 435-366mm Lengths

### Specifications

Cat. No.	Capacity (mL)	Graduation Interval (mL)	Hole Diameter x Length (mm)	Color Band	Pack Qty.	Case Qty.
WideTip, Two Lengths						
20200-11A	1	0.2	1.5 x 350	Red	10	90
20200-31A	3	0.5	2.0 x 350	Green	10	90
20200-51A	5	0.5	2.0 x 350	Blue	10	60
20200-101A	10	1	3.0 x 350	Orange	10	120
20200-12A	1	0.2	1.5 x 205	Red	10	90
20200-32A	3	0.5	2.0 x 205	Green	10	90
20200-52A	5	0.5	2.0 x 205	Blue	10	60
20200-102A	10	1	3.0 x 205	Orange	10	120
Large Hole, Short Length						
20210-52A	5	0.5	3.0 X 205	Blue	10	60
20210-102A	10	0.2	3.5 X 205	Orange	10	120
Graduated to the Tip						
20140-1A	1	0.01	1.2 X 360	N/A	10	60
20140-2A	2	0.02	1.2 X 360	N/A	10	60
20140-5A	5	0.05	1.2 X 366	N/A	10	60
20140-10A	10	0.1	1.4 X 366	N/A	10	60
20140-20A	20	0.1	1.5 X 435	N/A	10	60



Available  
at  
Comparable  
Price



**WARNING :** This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to [www.P65warnings.ca.gov](http://www.P65warnings.ca.gov)

**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER

1-800-DAWN-SCI

FOR INQUIRY

+1 732-902-6300

SALES@DAWNSCIENTIFIC.COM