



DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



Measuring Pipets

Class A, Reusable Glass, Color-Coded, w/ Colored Markings, To Deliver

SIBATA 2010A-Series Measuring Pipets meet ASTM E-1293 Style 1, Class A Specifications & USP Standards for Volumetric Glassware.

These pipets are graduated to the baseline above the tip (tempered) and are also calibrated to deliver. Feature permanent amber markings and color-codings which are enameled onto the glass.

The top end of sizes 5mL to 25mL is constricted.

Specifications

Cat. No.	Capacity (mL)	Tolerance (\pm mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Color Code
2010A-01-c	0.1	0.0025	0.01	6x300	White
2010A-02-c	0.2	0.004	0.01	6x300	Black
2010A-1-d	1.0	0.01	0.1	7x360	Red
2010A-1-c	1.0	0.01	0.01	7x360	Yellow
2010A-2-d	2.0	0.01	0.1	7x360	Green
2010A-5	5.0	0.02	0.1	8x365	Blue
2010A-10	10.0	0.03	0.1	10x370	Orange
2010A-25	25.0	0.05	0.1	13x455	White

Packaging

Cat. No.	Pack Qty.	Case Qty.
2010A-01-c	12	72
2010A-02-c	12	72
2010A-1-d	12	72
2010A-1-c	12	72
2010A-2-d	12	72
2010A-5	12	72
2010A-10	12	72
2010A-25	12	72

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

WARNING : This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to www.P65warnings.ca.gov

Available
at
Compatible
Price



DAWN SCIENTIFIC INC

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER
1-800-DAWN-SCI

FOR INQUIRY
+1 732-902-6300
SALES@DAWNSCIENTIFIC.COM

