



**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER



# Serological Pipets

Class A, Reusable Glass, Color-Coded, w/ Colored Markings, To Deliver

SIBATA 2011A-Series Serological Pipets meet ASTM E-1044 Style 1 Specifications.

These pipets are graduated to the tip (tempered) and are also calibrated to deliver. Feature permanent amber markings and color-codings which are enameled onto the glass.

The top end of sizes 5mL to 25mL is constricted.

## Specifications

Cat. No.	Capacity (mL)	Tolerance (±mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Color Code
2011A-01-c	0.1	0.0025	0.01	6 x 300	White
2011A-02-c	0.2	0.004	0.01	6 x 300	Black
2011A-1-d	1.0	0.01	0.1	7 x 360	Red
2011A-1-c	1.0	0.01	0.01	7 x 360	Yellow
2011A-2-d	2.0	0.01	0.1	7 x 360	Green
2011A-5	5.0	0.02	0.1	8 x 365	Blue
2011A-10	10.0	0.03	0.1	10 x 370	Orange
2011A-25	25.0	0.05	0.1	13 x 455	White

## Packaging

Cat. No.	Pack Qty.	Case Qty.
2011A-01-c	12	72
2011A-02-c	12	72
2011A-1-d	12	72
2011A-1-c	12	72
2011A-2-d	12	72
2011A-5	12	72
2011A-10	12	72
2011A-25	12	72



Available at Comparable Price

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

**WARNING :** This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to [www.P65warnings.ca.gov](http://www.P65warnings.ca.gov)



**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER  
**1-800-DAWN-SCI**

FOR INQUIRY  
**+1 732-902-6300**  
SALES@DAWNSCIENTIFIC.COM