



## Volumetric Pipets

**Class A, Reusable Glass, Color-Coded, w/ Colored Markings, To Deliver**

SIBATA 2040A-Series Volumetric Pipets meet ASTM E-969 Class A Specifications & USP Standards for Volumetric Glassware.

These pipets have sharp and permanent graduation lines and are calibrated to deliver. Feature permanent amber markings and color-codings which are enameled onto the glass.

### Specifications

Cat. No.	Capacity (mL)	Tolerance (±mL)	Approx. Length (mm)	Color Code	Pack Qty.	Case Qty.
2040A-05	0.5	0.005	300	2-Black	12	72
2040A-1	1.0	0.006	345	Blue	12	72
2040A-2	2.0	0.006	355	Orange	12	72
2040A-3	3.0	0.01	360	Black	12	72
2040A-4	4.0	0.01	365	2-Red	12	72
2040A-5	5.0	0.01	375	White	12	72
2040A-6	6.0	0.01	390	2-Orange	12	72
2040A-7	7.0	0.01	410	2-Green	12	72
2040A-8	8.0	0.02	430	Blue	12	72
2040A-9	9.0	0.02	436	Black	12	72
2040A-10	10.0	0.02	440	Red	12	72
2040A-13	13.0	0.03	460	White	6	36
2040A-15	15.0	0.03	460	Green	6	36
2040A-20	20.0	0.03	475	Yellow	6	36
2040A-25	25.0	0.03	480	Blue	6	36
2040A-50	50.0	0.05	540	Red	6	36
2040A-100	100.0	0.08	600	Yellow	6	36

"Dawn Scientific Inc is women owned small business located in central New Jersey. We have been serving our esteemed customers from various industries since 1978. Our priority is to supply quality lab products and chemicals from top of the range suppliers to our customers in Pharmaceuticals, Biopharma, Chemical, Cosmetics, Dietary & Herbal companies and many more whilst keeping in mind that better prices need not be compromised. We take great pride in serving other small businesses as well as educational institutes."

**WARNING :** This product can expose you to chemicals including lead, which is known to the state of California to cause cancer & birth defects or other reproductive harm. For more information go to [www.P65warnings.ca.gov](http://www.P65warnings.ca.gov)

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER



Available  
at  
Comparable  
Price



**DAWN SCIENTIFIC INC**

121 LIBERTY STREET, METUCHEN, NJ, 08840

TOLL-FREE NUMBER  
**1-800-DAWN-SCI**

FOR INQUIRY  
**+1 732-902-6300**  
SALES@DAWNSCIENTIFIC.COM