

**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

New

Product Announcement



BioWand

Personal UV sanitizer



PCR Plates



PPE/Masks



Micro Tubes



Dissection Tools

UV-C Light

Destroys
up to**99%**

of viruses and bacteria

New!

WWW.DAWNSCIENTIFIC.COM

Cat-No.	Description
U1001-W	Bio-Wand handheld personal sanitizer

Product Features

- Portable UV Disinfection
- Kills upto 99% of bacteria and viruses
- Destroys DNA & RNA
- Handheld, cordless, rechargeable

MTC's new BioWand uses UV light (265 to 280nm) to disinfect any object/surface in the lab and beyond. Bio-Wand can be set for a standard one-minute cleaning cycle or a more intense three-minute cycle. Simple press the power/ start button and hold the light 1-2 inches from the surface and slowly move it over the entire target area. Allow the UV light to shine on the desired area for 5 to 10 seconds to ensure optimum exposure. The Bio-Wand provides 2 hours of perating time and includes a USB-C cable for charging and portability.

BioWand Personal UV Sanitizer

UV Light Source	: LED, UVC 265-280nm
Power Input	: USB 5Volt, 2Amp
Rated Power	: 2W
Battery	: Lithium 3.7v, 500mAh
Operating Time	: 2 Hours per charge
Cycle Times	: 60 Sec, 180 sec

WWW.DAWNSCIENTIFIC.COM