



DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER



Exclusively
available at
Dawn Scientific INC

CUVETTE

Part No. DS-29500-004

✓ Material

Customized Large Glass Cuvette
It is exclusively made up of Borosilicate glass.

✓ Uses

Specially designed to hold the samples for the spectroscopic analysis.

✓ Construction

It has all sides fused, No Adhesives.
Dimensions are 54mm X 54mm, with 2mm thickness

✓ Features

Path length - 10mm
Path length Tolerance - ± 0.01 mm