




**DAWN SCIENTIFIC INC**

YOUR SCIENTIFIC PARTNER

**NEW  
PRODUCT**



**Part No. DS-91000-222**

# ELECTROMAGNETIC COIL

> Provides horizontal movement of the swinging funnels in automatic distillation heads for reflux or take off.

**Exclusively  
Available at  
DAWN SCIENTIFIC INC**

- > Can be operated for twelve weeks straight without letup and without failures. There is no burn-out and no significant heat build up to ignition.
- > Specially designed for use with automatic times (optional). This coil is a potted type horseshoe core design to operate large automatic distilling heads, up to 3 1/4" (83mm) O.D. with adjustable clamp.