

DAWN SCIENTIFIC INC

YOUR SCIENTIFIC PARTNER

HPLC Solvents

Gradient Grade Solvents for your critical HPLC analysis



Product Features

- Performance are validated by Gradient Elution test at 210 nm & 254 nm for Acetonitrile & 230 nm & 254 nm for Methanol
- Our Gradient elution method complies to Reagent Chemicals eleventh addition.
- Our Gradient solvents are in the scope of ISO17025
- Elution test on specified wavelength can be conducted on request

Leading Manufactures

MILLIPORE SIGMA

Alfa Aesar
by Thermo Fisher Scientific

TGI

GREENFIELD
GLOBAL
LABORATORY SUPPLIES & EQUIPMENT

CSI

M TEDIA
High Purity Solvents

45
CENTRE

Product Range

- » Acetic acid glacial HPLC & Spectroscopy
- » Acetone for HPLC & Spectroscopy
- » Acetonitrile For HPLC & Spectroscopy
- » Acetonitrile Gradient Grade for HPLC
- » Benzene for HPLC & Spectroscopy
- » Benzyl alcohol for HPLC & Spectroscopy
- » Butan-1 ol for HPLC & Spectroscopy
- » tert-Butanol for HPLC & Spectroscopy
- » Butan-2-one for HPLC & Spectroscopy
- » tert-Butyl methyl ether for HPLC
- » Chloroform HPLC & Spectroscopy
- » Cyclohexane HPLC & Spectroscopy
- » 1,2-Dichloroethane HPLC & Spectroscopy
- » Dichloromethane HPLC & Spectroscopy
- » Diethylamine for HPLC & Spectroscopy
- » Diethyl ether for HPLC & Spectroscopy
- » 1,2 Dimethoxyethane for HPLC & Spectroscopy
- » N,N-Dimethyl acetamide for HPLC & Spectroscopy
- » N,N-Dimethyl formamide for HPLC & Spectroscopy
- » 1,4-Dioxan for HPLC & Spectroscopy
- » Dimethyl sulphoxide for HPLC & Spectroscopy



DAWN SCIENTIFIC INC

FOR INQUIRY :+1 732-902-6300
sales@dawnscientific.com

TOLL-FREE NO.1-800-DAWN-SCI
www.dawnscientific.com