

Product : Benzoic acid
Catalogue : CH2537
CAS : 65-85-0
Formula : C₆H₅COOH
MW : 122.12
Grade : ACS

Test	Specification
Appearance	Colorless or white solid
Assay	≥ 99.5%
Freezing point	122.0°C to 123.0°C
Residue after ignition	≤ 0.005%
Insoluble in methanol	≤ 0.005%
Chlorine compounds (as Cl)	≤ 0.005%
Sulfur compounds (as S)	≤ 0.002%
Heavy metals (as Pb)	≤ 5 ppm
Substance reducing permanganate	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

