

Product : Ammonium Carbonate, ACS
Catalog# : CH1030
CAS : 8000-73-5
Formula : $(\text{NH}_4)_2\text{CO}_3$
MW : 96.06
Grade : ACS

TEST	SPECIFICATION
Appearance	White powder with odor of ammonia
Assay	$\geq 30.0 \%$
Insoluble matter	$\leq 0.005 \%$
Nonvolatile matter	$\leq 0.01 \%$
Chloride (Cl)	$\leq 5 \text{ ppm}$
Sulfur compounds (as SO ₄)	$\leq 0.002 \%$
Heavy metals	$\leq 5 \text{ ppm}$
Iron (Fe)	$\leq 5 \text{ ppm}$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

