

Product : Ammonium Chloride Lab
Catalog# : CH1032
CAS : 12125-02-9
Formula : NH₄Cl
MW : 53.49
Grade : Lab

TEST	SPECIFICATION
Appearance	White solid
Assay	≥ 99.4 %
pH of 5 % solution at 25°C	4.5 – 5.5
Insoluble matter	≤ 0.05 %
Residue after ignition	≤ 0.01 %
Calcium (Ca)	≤ 0.01 %
Magnesium (Mg)	≤ 5 ppm
Phosphate (PO ₄)	≤ 2 ppm
Sulfate (SO ₄)	≤ 0.02 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

