

Product : Ammonium Fluoride ACS
Catalog# : CH1038
CAS : 12125-01-8
Formula : NH₄F
MW : 37.04
Grade : ACS

TEST	SPECIFICATION
Appearance	Colorless or white, deliquescent solid
Assay	≥ 98.0 %
Insoluble matter	≤ 0.005 %
Residue after ignition	≤ 0.01 %
Chloride (Cl)	≤ 0.001 %
Sulfate (SO ₄)	≤ 0.005 %
Heavy metals (as Pb)	≤ 5 ppm
Iron (Fe)	≤ 5 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

