

Product : Ammonium Thiocyanate ACS
Catalog# : CH1055
CAS : 1762-95-4
Formula : NH₄SCN
MW : 76.12
Grade : ACS

TEST	SPECIFICATION
Appearance	Colorless, deliquescent liquid
Assay	≥ 97.5 %
pH of a 5 % solution at 25.0°C	4.5 – 6.0
Insoluble matter	≤ 0.005 %
Residue after ignition	≤ 0.025 %
Chloride (Cl)	≤ 0.005 %
Sulfate (SO ₄)	≤ 0.005 %
Heavy metals	≤ 5ppm
Iron (Fe)	≤ 3 ppm
Iodine-consuming substance	≤ 0.004 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

