

Product : Arsenic Trioxide ACS  
Catalog# : CH1060  
CAS : 1327-53-3  
Formula : As<sub>2</sub>O<sub>3</sub>  
MW : 197.84  
Grade : ACS

TEST	SPECIFICATION
Appearance	White solid
Assay	99.95 – 100.05 %
Residue after ignition	≤ 0.02 %
Insoluble in dilute hydrochloric acid	≤ 0.01 %
Chloride (Cl)	≤ 0.005 %
Sulfide (S)	Passes test
Antimony (Sb)	≤ 0.05 %
Lead (Pb)	≤ 0.002 %
Iron (Fe)	≤ 5 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

