

Product : Boric Acid Granular ACS
Catalog# : CH1090
CAS : 10043-35-3
Formula : H₃BO₃
MW : 61.83
Grade : ACS

TEST	SPECIFICATION
Appearance	Colorless or white solid
Assay	≥ 99.5 %
Insoluble in methanol	≤ 0.005 %
Nonvolatile with methanol	≤ 0.05 %
Chloride (Cl)	≤ 0.001 %
Phosphate (PO ₄)	≤ 0.001 %
Sulfate (SO ₄)	≤ 0.01 %
Heavy metals (as Pb)	≤ 0.001 %
Iron (Fe)	≤ 0.001 %
Calcium (Ca)	≤ 0.005 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

