

Product : Boric Acid Lab
Catalog# : CH1091
CAS : 10043-35-3
Formula : H_3BO_3
MW : 61.83
Grade : Lab

TEST	SPECIFICATION
Appearance	Colorless or white solid
Assay	$\geq 99.5\%$
Insoluble in methanol	$\leq 0.005\%$
Nonvolatile with methanol	$\leq 0.05\%$
Chloride (Cl)	$\leq 0.001\%$
Phosphate (PO_4)	$\leq 0.001\%$
Sulfate (SO_4)	$\leq 0.01\%$
Calcium (Ca)	$\leq 0.005\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

