

Product : Ferrous Ammonium Sulfate Hexahydrate
 Catalog# : CH1232
 CAS : 7783-85-9
 Formula : $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$
 MW : 392.13
 Grade : Lab

TEST	SPECIFICATION
Appearance	Pale blue – green solid
Assay	98.5 – 101.5 %
Insoluble matter	≤ 0.01 %
Phosphate (PO ₄)	≤ 0.003 %
Calcium (Ca)	≤ 0.005 %
Copper (Cu)	≤ 0.003 %
Magnesium (Mg)	≤ 0.002 %
Manganese (Mn)	≤ 0.01 %
Potassium (K)	≤ 0.002 %
Sodium (Na)	≤ 0.02 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

