



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Alizarin Red S 1%
Catalog:	IS10753
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Alizarin Red S Monohydrate, Sodium Salt	130-22-3	Certified

Specification:

Test	Specification
Appearance	Dark red liquid
Suitability for Use	Yellow (pH 4.6) - Pink (pH 6.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

