



Product Specification

Product Information:

Products:	Bromocresol Green 0.04%
Catalog:	IS10754
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Bromocresol Green, Sodium Salt	62625-32-5	ACS

Specification:

Test	Specification
Appearance	Dark blue liquid
Suitability for Use	Blue (pH >5.4) - Green (pH 4.5) -Yellow (pH < 4.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

