



Product Specification

Product Information:

Products:	Bromocresol Green 0.1% (in IPA)
Catalog:	IS10755
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Isopropyl Alcohol	67-63-0	ACS
Bromocresol Green, sultone	76-60-8	ACS

Specification:

Test	Specification
Appearance	Orange liquid
Suitability for Use	Blue (pH 5.4) - Yellow (pH 3.8)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

