



# Product Specification

## Product Information:

<b>Products:</b>	Bromocresol Green/Methyl Red in IPA
<b>Catalog:</b>	IS10756
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Isopropyl Alcohol	67-63-0	ACS
Bromocresol Green, Sodium Salt	62625-32-5	ACS
Methyl Red, Sodium Salt	845-10-3	ACS

## Specification:

Test	Specification
Appearance	Dark colored liquid
Suitability for Use	Greenish-Blue (pH 5.4) – Light Pink (pH 4.5), with Gray warning color

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

