



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Bromophenol Blue 1%
Catalog:	IS10764
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Bromophenol Blue, Sodium Salt	62625-28-9	ACS

Specification:

Test	Specification
Appearance	Dark reddish-purple liquid
Suitability for Use	Yellow (pH < 3.0) - Green (pH 3.4to 3.7) - Blue (pH >4.6)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

