



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Bromothymol Blue 0.04%
Catalog:	IS10765
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Bromothymol Blue Sodium Salt	34722-90-2	ACS

Specification:

Test	Specification
Appearance	Dark colored liquid
pH at 25°C	6.7-6.9
Suitability for Use	Yellow (pH < 6.0) - Green (pH 6.7) - Blue (pH >7.6)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

