



Product Specification

Product Information:

| | |
|------------------|--------------------|
| Products: | Calmagite 0.05% |
| Catalog: | IS10769 |
| Grade: | Indicator solution |

Reagent Used:

| Reagents | CAS # | Grade |
|-----------|-----------|-----------------|
| Water | 7732-18-5 | ACS/ASTM/USP/EP |
| Calmagite | 3147-14-6 | Indicator |

Specification:

| Test | Specification |
|---------------------|---|
| Appearance | Red liquid |
| Suitability for Use | Wine Red (excess Ca/Mg) – Pure Blue (excess EDTA) |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

