



# Product Specification

## Product Information:

<b>Products:</b>	Eosin Y 1% in Ethanol
<b>Catalog:</b>	IS10776
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Ethyl Alcohol	64-17-5	Reagent
Eosin Y	17372-87-1	Certified

## Specification:

Test	Specification
Appearance	Murky green liquid
UV/Visible Spectrum	To Pass Test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

