



Product Specification

Product Information:

Products:	meta-Cresol Purple 0.1%
Catalog:	IS10785
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
m-Cresol Purple, Sodium Salt	62625-31-4	Reagent

Specification:

Test	Specification
Appearance	Dark, reddish-purple
Visual transition range	1.2 - 2.8, red to yellow (acid) 7.4 - 9.0, yellow to violet (alkaline)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

